

## Product datasheet for **RC234777**

### **HSP90AB1 (NM\_001271969) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	HSP90AB1 (NM_001271969) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	HSP90AB1
Synonyms:	D6S182; HSP84; HSP90B; HSPC2; HSPCB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC234777 representing NM\_001271969  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCCTGAGGAAGTGCACCATGGAGAGGAGGAGGTGGAGACTTTTGCCTTCAGGCAGAAATGCCCAAC  
 TCATGTCCCTCATCATCAATACCTTCTATTCCAACAAGGAGATTTTCTTCGGGAGTTGATCTCTAATGC  
 TTCTGATGCCTTGACAAGATTGCTATGAGAGCCTGACAGACCCCTCGAAGTTGGACAGTGGTAAAGAG  
 CTGAAAATTGACATCATCCCAACCCCTCAGGAACGTACCCTGACTTTGGTAGACACAGGCATTGGCATGA  
 CCAAAGCTGATCTCATAAATAATTTGGGAACCTTGCCAAGTCTGGTACTAAAGCATTATGGAGGCTCT  
 TCAGGCTGGTGCAGACATCTCCATGATTGGGCAGTTTGGTGTGGCTTTTATTCTGCCTACTTGGTGGCA  
 GAGAAAGTGGTTGTGATCACAAGCACAACGATGATGAACAGTATGCTTGGGAGTCTTCTGCTGGAGGTT  
 CCTTCACTGTGCGTGCTGACCATGGTGAAGCCATTGGCAGGGTACCAAAGTGATCTCCATCTTAAGA  
 AGATCAGACAGAGTACCTAGAAGAGAGGGCGGTCAAAGAAGTAGTGAAGAAGCATTCTCAGTTCATAGGC  
 TATCCCATCACCTTTATTTGGAGAAGGAACGAGAGAAGGAAATTAGTGATGATGAGGCAGAGGAAGAGA  
 AAGGTGAGAAAAGAGGAAGATAAAGATGATGAAGAAAAGCCCAAGATCGAAGATGTGGGTTTCAGATGA  
 GGAGGATGACAGCGGTAAGGATAAGAAGAAGAAAATAAGAAGATCAAAGAGAAATACATTGATCAGGAA  
 GAATAACAAGACCAAGCCTATTTGGACCAGAAACCCTGATGACATCACCAAGAGGAGTATGGAGAAT  
 TCTACAAGAGCCTCACTAATGACTGGGAAGACCACTTGGCAGTCAAGCACTTTTCTGTAGAAGGTGAGTT  
 GGAATTCAGGGCATTGCTATTTATCTCGTGGGCTCCCTTTGACCTTTTGGAGAACAAGAAGAAAAAG  
 AACAACATCAAACCTCTATGTCCCGGTGTGTTTCATCATGGACAGCTGTGATGAGTTGATACCAGAGTATC  
 TCAATTTTATCCGTGGTGTGGTTGACTCTGAGGATCTGCCCTGAACATCTCCGAGAAATGCTCCAGCA  
 GAGCAAAATCTTGAAAGTATTGCAAAAACATTTGTTAAGAAGTGCCTTGAGCTCTTCTGAGCTGGCA  
 GAAGACAAGGAGAATTACAAGAAATCTATGAGGCATTCTCAAAAATCTCAAGCTTGAATCCACGAAG  
 ACTCCACTAACCGCCCGCCTGTCTGAGCTGCTGCGCTATCATACCTCCAGTCTGGAGATGAGATGAC  
 ATCTCTGTGAGATATGTTTCTCGCATGAAGGAGACACAGAAGTCCATCTATTACATCACTGGTGAAGC  
 AAAGAGCAGTGGCCAACCTCAGCTTTTGTGGAGCGAGTGCAGAAACGGGCTTCGAGGTGGTATATGA  
 CCGAGCCATTGACGAGTACTGTGTGCAGCAGCTCAAGGAATTTGATGGGAAGAGCCTGGTCTCAGTTAC  
 CAAGGAGGCTGGAGCTGCCTGAGGATGAGGAGGAGAAGAAGAAGATGGAAGAGAGCAAGGCAAAGTTT  
 GAGAACCTCTGCAAGCTCATGAAAGAACTTAGATAAGAAGGTTGAGAAGGTGACAATCTCCAATAGAC  
 TTGTGTCTTACCTTGCTGCATTGTGACCAGCACCTACGGCTGGACAGCCAATATGGAGCGGATCATGAA  
 AGCCAGGCATTCGGGACAACCTCACCATGGGCTATATGATGGCCAAAAGCACCTGGAGATCAACCCCT  
 GACCACCCATTGTGGAGACGCTGCGGCAGAAGGCTGAGGCCGACAAGAATGATAAGGCAGTTAAGGACC  
 TGGTGGTGTGCTGTTTGAACCGCCCTGCTATCTTCTGGCTTTTCCCTTGAGGATCCCCAGACCCACTC  
 CAACCGCATCTATCGCATGATCAAGCTAGGTCTAGGTATTGATGAAGATGAAGTGGCAGCAGAGGAACCC  
 AATGCTGCAGTTCCTGATGAGATCCCCCTCTCGAGGGCGATGAGGATGCGTCTCGCATGGAAGAAGTCG  
 AT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC234777 representing NM\_001271969  
Red=Cloning site Green=Tags(s)


```
MPEEVHHGEEEVETFAFQAEIAQLMSLIINTFYNSKEIFLRELISNASDALDKIRYESLTDPSKLDSGKE
LKIDIIIPNPQERTLTLVDTGIGMTKADLINNLTIAKSGTKAFMEALQAGADISMGQFVGVFYSAYLVA
EKVVVITKHNDEQYAWESSAGGSFTVRADHGPEIGRGTIVILHLKEDQTEYLEERRVKEVVKHSQFIG
YPITLYLEKEREKEISSDDEAEKEEKEEEDKDEEKPKIEDVGSDEEDDSGDKKKKTKIKEKYIDQE
ELNKTPIWTRNPDDITQEEYGEFYKSLTNDWEDHLAVKHFVSEGQLEFRALLFIPRRAPFDLFENKKKK
NNIKLYVRRVIMDSCDELIPEYLNFRGVVDSDELPLNISREMLQQSKILKVIKNIIVKCCLELFSELA
EDKENYKFFYEAFSKNLKLGIHEDSTNRRRLSELLRYHTSQSGDEMTSLSEYVSRMKETQKSIYYITGES
KEQVANSFAVERVRKRGFEVVMTEPIDEYCVQQLKEFDGKSLVSVTKEGLELPEDEEEKKKMEESKAKF
ENLCKLMKEILDKKVEKVTISNRLVSSPCIVTSTYGTWANMERIMKAQALRDNSTMGYMAKKHLEINP
DHPIVETLRQKAEADKNDKAVKDLVLLFETALLSSGFSLEDQPQTHSNRIYRMIKLGLGIDEVEVAEEP
NAAVPDEIPPLEGDEDASRMEEVD
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mg2549\\_g07.zip](https://cdn.origene.com/chromatograms/mg2549_g07.zip)

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001271969


**ORF Size:** 2172 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001271969.1</a></u> , <u><a href="#">NP_001258898.1</a></u>
<b>RefSeq Size:</b>	2716 bp
<b>RefSeq ORF:</b>	2175 bp
<b>Locus ID:</b>	3326
<b>UniProt ID:</b>	<u><a href="#">P08238</a></u>
<b>Cytogenetics:</b>	6p21.1
<b>Protein Families:</b>	Druggable Genome, Stem cell - Pluripotency
<b>Protein Pathways:</b>	Antigen processing and presentation, NOD-like receptor signaling pathway, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer
<b>MW:</b>	83.3 kDa
<b>Gene Summary:</b>	This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes. [provided by RefSeq, Dec 2012]

Product images:



Circular map for RC234777