

Product datasheet for **RC234797**

BRSK2 (NM_001256627) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRSK2 (NM_001256627) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BRSK2
Synonyms:	C11orf7; PEN11B; SAD1; SADA; STK29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234797 representing NM_001256627
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACATCGACGGGAAGGACGGCGCGCGCAGCACGCGCAGTATGTTGGGCCCTACCGGTGGAGAAGA
 CGCTGGGCAAGGGCAGACAGGTCTGGTGAAGCTGGGGTTCACTGCGTCACCTGCCAGAAGGTGGCCAT
 CAAGATCGTCAACCGTGAGAAGCTCAGCGAGTCCGGTCTGATGAAGGTGGAGCGGGAGATCGCGATCCTG
 AAGCTCATTGAGCACCCACGTCCTAAAGCTGCACGACGTTTATGAAAACAAAAATATTTGTACCTGG
 TGCTAGAACACGTGTCAGGTGGTGAAGCTCTCGACTACCTGGTGAAGAAGGGGAGGCTGACGCCTAAGGA
 GGCTCGGAAGTTCTCCGGCAGATCATCTCTGCGCTGGACTTCTGCCACAGCCACTCCATATGCCACAGG
 GATCTGAAACCTGAAACCTCCTGCTGGACGAGAACAACATCCGCATCGCAGACTTTGGCATGGCGT
 CCCTGCAGTTGGCGACAGCCTGTTGGAGACCAGCTGTGGTCCCCCACTACGCCTGCCCGAGGTGAT
 CCGGGGGGAGAAGTATGACGGCCGAAGGCGGACGTGTGGAGCTCGGGCGTCATCTGTTCGCTTGCTG
 GTGGGGCTCTGCCCTTCGACGATGACAACCTTGCAGACGCTGCTGGAGAAGGTGAAGCGGGCGTGTTC
 ACATGCCGCACCTTATCCCGCCGACTGCCAGAGTCTGCTACGGGGCATGATCGAGGTGGACGCCGACG
 CCGCTCACGCTAGAGCACATTCAGAAACATATGGTATATAGGGGGCAAGAATGAGCCGAACCAGAG
 CAGCCATTCTCGAAGGTGCAGATCCGCTCGCTGCCAGCCTGGAGGACATCGACCCGACGTGCTGG
 ACAGCATGCACTCACTGGGCTGCTCCGAGACCGCAACAAGTCTGTCAGGACCTGCTGTCCGAGGAGGA
 GAACCAGGAGAAGATGATTTACTTCTCCTCTGGACCGAAAGAAAGGTACCCGAGCCAGGAGGATGAG
 GACCTGCCCCCGGAACGAGATAGACCTCCCCGGAAGCGTGTGGACTCCCCGATGCTGAACCGGCACG
 GCAAGCGCGGCCAGAACGCAAATCCATGGAGGTGCTCAGCGTGACGGACGGCGGCTCCCCGGTGCCTGC
 GCGGCGGCCATTGAGATGGCCACGACGCGCCAGAGGTCTCGGTCCATCAGCGGTGCCTCCTCAGGCCTT
 TCCACCAGCCACTCAGCAGCCCCGGGTACCCCTACCCCTCACCAAGGGGCAGTCCCCCTCCCCACCC
 CCAAGGGGACACCTGTCCACAGCCAAAGGAGAGCCCGGCTGGCAGCCCAACCCACGCCCCCGTCCAG
 CCCCAGCGTCGGAGGGGTGCCCTGGAGGGCGGGCTCAACTCCATCAAGAACAGCTTTCTGGGCTCACCC
 CGCTTCCACCGCCGAAACTGCAAGTCCGACGCGGAGGAGATGTCAAACCTGACACCAGAGTCGTCCC
 CAGAGCTGGCGAAGAAGTCTGGTTTGGAACTTATCAGCCTGGAGAAGGAGGAGCAGATCTTCGTGGT
 CATCAAAGACAACTCTGAGCTCCATCAAGGCTGACATCGTGCACGCCTTCTGTGATCCCAGTCTC
 AGCCACAGCGTCACTCCCAAACGAGCTTCCGGGCGGAGTACAAGGCCACGGGGGGCCAGCCGTGTTCC
 AGAAGCCGGTCAAGTCCAGGTTGATATCACTACACGGAGGGTGGGAGGCGCAGAAGGAGAACCGCAT
 CTACTCCGTACCTTACCCCTGCTCTCAGGCCACGCCGTCGCTTCAAGAGGGTGGTGGAGACCATCCAG
 GCCAGCTGCTGAGCACACAGCCCGCTGCGGCCAGCACTTGTGAGACACCACTAAGTATGGAAA
 TGATGACGGGGCGGCTTTCCAAATGTGGCAGCCATTGAGTAACCTTTTGGACGTAATTAACAACCTTTT
 TTCAGACGAGAAGAACGGGCAGGCGGCCAGGCCCCAGCACGCCCGCAAGCGGAGTGCCACGGCCCA
 CTCGGTGACTCCGCGCCGCTGGCCCTGGCCCCGAGGGGACGCCGAGTACCCAACGGGCAAGGACACGG
 CCAAGATGGGCCCGCCACCGCCCGCGGAGCAGCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC234797 representing NM_001256627
Red=Cloning site Green=Tags(s)

MTSTGKDGGAQHAQYVGPYRLEKTLGKGQTGLVKLVHCVTCQKVAIKIVNREKLSSEVLMKVEREIAIL
 KLIEHPHVLKLDVYENKKYLVLVLEHVS GGELFDYL VKKGR LTPKEARKFFRQIISALDFCHSHSICHR
 DLKPENLLLDEKNNIRIADFGMASLQVGD S LLETSCGSPHYACPEVIRGEKYDGRKADVWSCGVILFALL
 VGALPFDDDLRQLLEKVKRGVFMHPHFIPDDCQSLLRGMI EVD AARRL TLEHIQKHIWYIGGKNEPEPE
 QPIPRKVQIRSLPSLEDIDPDVLD SMHSLGCFRDRNKLLQDLLSEENQEKM IYFLLLDRKERYPSQEDE
 DLPPRNEIDPPRKRVDSPMLNRHGKRRRPERKSMEVLSVTDGGSPVPARRA IEMAQHGRSRISGASSGL
 STSPLSSPRVTPHPSPRGSPLPTPKGTPVHTPKESPAGTPNPTPPSSPSVGGVPWRARLNSIKNSFLGSP
 RFHRRKLQVPTPEEMSNLTPESP ELAKKSWFGNFISLEKEEQIFVVIKDKPLSSIKADIVHAF LSIPLS
 SHSVISQTSFRAEYKATGGPAVFQKPVKFQVDITYTEGGEAQKENG IYSVTFLLSGPSRRFRKRVETIQ
 AQLLSTHDPAAQHLSDTTNCMEMMTGRLSKCGSPLSNFFDVIKQLFSDEKNGQAAQAPSTPAKRS AHGP
 LGDSAAAGPGGGDAEYPTGKDTAKMGPP TARREQP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001256627

ORF Size: 2208 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256627.2
RefSeq Size:	4563 bp
RefSeq ORF:	2211 bp
Locus ID:	9024
UniProt ID:	Q8IWQ3
Cytogenetics:	11p15.5
Protein Families:	Druggable Genome, Protein Kinase
MW:	82.1 kDa
Gene Summary:	<p>Serine/threonine-protein kinase that plays a key role in polarization of neurons and axonogenesis, cell cycle progress and insulin secretion. Phosphorylates CDK16, CDC25C, MAPT/TAU, PAK1 and WEE1. Following phosphorylation and activation by STK11/LKB1, acts as a key regulator of polarization of cortical neurons, probably by mediating phosphorylation of microtubule-associated proteins such as MAPT/TAU at 'Thr-529' and 'Ser-579'. Also regulates neuron polarization by mediating phosphorylation of WEE1 at 'Ser-642' in postmitotic neurons, leading to down-regulate WEE1 activity in polarized neurons. Plays a role in the regulation of the mitotic cell cycle progress and the onset of mitosis. Plays a role in the regulation of insulin secretion in response to elevated glucose levels, probably via phosphorylation of CDK16 and PAK1. While BRSK2 phosphorylated at Thr-174 can inhibit insulin secretion (PubMed:22798068), BRSK2 phosphorylated at Thr-260 can promote insulin secretion (PubMed:22669945). Regulates reorganization of the actin cytoskeleton. May play a role in the apoptotic response triggered by endoplasmic reticulum (ER) stress.</p> <p>[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for RC234797