

Product datasheet for **RC234798**

EZH2 (NM_001203248) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EZH2 (NM_001203248) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	EZH2
Synonyms:	ENX-1; ENX1; EZH2b; KMT6; KMT6A; WVS; WVS2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234798 representing NM_001203248
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCCAGACTGGGAAGAAATCTGAGAAGGGACCAGTTTGTGGCGGAAGCGTGTAAAATCAGAGTACA
 TGCAGACTGAGACAGCTCAAGAGTTTCAGACGAGCTGATGAAGTAAAGAGTATGTTTAGTTCCAATCGTCA
 GAAAAATTTGGAAAGAACGGAAATCTTAAACCAAGAATGGAAACAGCGAAGGATACAGCCTGTGCACATC
 CTGACTTCTTGTTCGGTGACCAGTGACTTGGATTTTCCAACACAAGTCATCCCATTAAGACTCTGAATG
 CAGTTGCTTCAGTACCCATAATGTATTCTTGGTCTCCCTACAGCAGAATTTTATGGTGGAAAGTAAAC
 TGTTTTACATAACATTCCTTATATGGGAGATGAAGTTTTAGATCAGGATGGTACTTTTCATTGAAGAACTA
 AAAAAAATTATGATGGGAAAGTACACGGGGATAGAGAATGTGGGTTTATAATGATGAAATTTTTGTGG
 AGTTGGTGAATGCCCTTGGTCAATATAATGATGATGACGATGATGATGATGGAGACGATCCTGAAGAAAG
 AGAAGAAAAGCAGAAAGATCTGGAGGATCACCGAGATGATAAAGAAAGCCGCCACCTCGGAAATTTCT
 TCTGATAAAATTTTTGAAGCCATTTCTCAATGTTCCAGATAAGGGCACAGCAGAAGAACTAAAGGAAA
 AATATAAAGAACTCACCGAACAGCAGCTCCAGGGCAGCTTCTCCTGAATGTACCCCCAACATAGATGG
 ACCAAATGCTAAATCTGTTTCAGAGAGAGCAAAGCTTACACTCCTTTCATACGCTTTTCTGTAGGCGATGT
 TTTAAATAGACTGCTTCTACATCCTTTTCATGCAACACCCAACACTTATAAGCGGAAGAACACAGAAA
 CAGCTCTAGACAACAACCTTGTGGACCACAGTGTACCAGCATTTGGAGGGAGCAAAGGAGTTTGTCTGC
 TGCTCTCACCGCTGAGCGGATAAAGACCCCAACAAACGTCAGGAGGCCGAGAAGAGGACGGCTTCCC
 AATAACAGTAGCAGGCCAGCACCCCAACATTAATGTGCTGGAATCAAAGGATACAGACAGTGAAGGG
 AAGCAGGGACTGAAACGGGGGGAGAGAACAATGATAAAGAAGAAGAAGAGAAGAAAGTAAACTTCGAG
 CTCTCTGAAGCAAATTCGCGGTGCAAAACACCAATAAAGATGAAGCAAATATTGAACCTCTGAGAAT
 GTGGAGTGGAGTGGTGTGAAGCCTCAATGTTTAGAGTCCCTATTGGCACTTACTATGACAATTTCTGTG
 CCATTGCTAGGTTAATTGGGACCAAAACATGTAGACAGGTGTATGAGTTTAGAGTCAAAGAATCTAGCAT
 CATAGCTCCAGCTCCCGTGAGGATGTGGATACTCCTCAAGGAAAAAGAAGAGGAAACACCGGTTGTGG
 GCTGCACACTGCAGAAAGATACAGCTGAAAAAGGACGGCTCCTTAACCATGTTTACAATATCAACCT
 GTGATCATCCACGGCAGCCTTGTGACAGTTCGTGCCCTTGTGTGATAGCACAAAATTTTTGTGAAAAGTT
 TTGTCAATGTAGTTCAGAGTGTCAAACCGCTTCCGGGATGCCGCTGCAAAGCACAGTGAACACCAAG
 CAGTGCCCGTGTACCTGGCTGTCCGAGAGTGTACCCTGACCTCTGTCTTACTTGTGGAGCCGCTGACC
 ATTGGGACAGTAAAAATGTGCTCTGCAAGAACTGCAGTATTCAGCGGGGCTCCAAAAAGCATCTATTGCT
 GGCACCATCTGACGTGGCAGGCTGGGGATTTTTATCAAAGATCCTGTGCAGAAAAATGAATTCATCTCA
 GAATACTGTGGAGAGATTATTTCTCAAGATGAAGCTGACAGAAGAGGGAAAGTGTATGATAAATACATGT
 GCAGCTTTCTGTTCAACTGAACAATGATTTTGTGGTGGATGCAACCCGCAAGGGTAACAAAATTCGTTT
 TGCAAATCATTCCGTAATCCAACTGCTATGCAAAAGTTATGATGGTTAACGGTGTACACAGGATAGGT
 ATTTTTGCCAAGAGGCCATCCAGACTGGCGAAGAGCTGTTTTTTGATTACAGATACAGCCAGGCTGATG
 CCCTGAAGTATGTCGGCATCGAAAGAGAAATGGAATCCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234798 representing NM_001203248
Red=Cloning site Green=Tags(s)


MGQTGKKSEKGPVCWRKRVKSEYMRLRQLKRFRADEVKSMFSSNRQKILERTEILNQEWKQRRIQPVHI
 LTSCSVTSDLDFPTQVIPLKTLNAVASVPIMYSWSPLQQNFMVEDETVLHNIPYMGDEVLDQDGTFFI EEL
 IKNYDGKLVHGDRECGFINDEIFVELVNALGQYNDDDDDDDDPEEREKQKDL EDRDDKESRPPRKF
 SDKIFEAISSMFPDKGTAEELKEKYKELTEQQLPGALPPECTPNIDGPNKSVQREQLHSFHTLFCRRC
 FKDYDCFLHPFHATPNTYKRKNTETALDNKPCGPQCYQHLEGAKEFAAALTAERIKTPPKRPGRRRGRLP
 NNSSRPSTPTINVLESKDTSDREAGTETGGENNDKEEEEKDETS SSSSEANSRCQTPIKMKPNIEPPEN
 VEWSGAEASMFVRLIGTYDNFCAIARLIGTKTRQVYEFVYKESIIAPAPAEDVTPPRKKRKHRLW
 AAHCRKIQLKKGSSNHVYNYQPCDHPRQPCDSSCPVIAQNFCEKFCQCSSECQNRFPGCCKAQCNTK
 QCPCYLAVRECDPDLCLTCGAADHWDSKNVSKNCISIQRGSKHLLLAPSDVAGWGIKIDPVQKNEFIS
 EYCGEII SQDEADRRGKVDKYMCSFLFNLNDFVVDATRKGNKIRFANHSVNPNCYAKVMMVNGDHRIG
 IFAKRAIQTGEELFFDYRYSQADALKYVGIEREMEIP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001203248

ORF Size: 2211 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001203248.2](#)

RefSeq Size: 2681 bp

RefSeq ORF: 2214 bp

Locus ID: 2146

UniProt ID: [Q15910](#)


Cytogenetics: 7q36.1

Protein Families: Druggable Genome, Transcription Factors

MW: 84.8 kDa

Gene Summary: This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Feb 2011]

Product images:



Circular map for RC234798