

Product datasheet for **RC234818**

TAX1BP1 (NM_001206901) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAX1BP1 (NM_001206901) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TAX1BP1
Synonyms:	CALCOCO3; T6BP; TXBP151
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC234818 representing NM_001206901
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGACATCCTTTCAAGAAGTCCCATTGCAGACTTCCAACCTTGCCCATGTCATCTTTCAAATGTGGCCA
 AGAGTTACCTTCCTAATGCACACCTGGAATGTCATTACACCTTAACTCCATATATTCATCCACATCCAAA
 AGATTGGGTTGGTATATCAAGGTTGGATGGAGTACTGCTCGTGATTATTACACGTTTTTATGGTCCCCT
 ATGCCTGAACATTATGTGGAAGGATCAACAGTCAATTGTGTACTAGCATTCCAAGGATATTACCTTCCAA
 ATGATGATGGAGAATTTATCAGTTCTGTTACGTTACCCATAAGGGTGAAATTCGTGGAGCAAGTACACC
 TTTCCAGTTTCGAGCTTCTTCTCCAGTTGAAGAGCTGCTTACTATGGAAGATGAAGGAAATTCGTGACATG
 TTAGTGGTGACCACAAAAGCAGGCTTCTTGAGTTGAAAATTGAGAAAACCATGAAAGAAAAGAAGAAC
 TGTTAAAGTTAATTGCCGTTCTGAAAAAGAAACAGCACAACTTCGAGAACAAGTTGGGAGAATGAAAAG
 AGAAGTTAACCATGAGAAAGAAAGATGTGACCAACTGCAAGCAGAACAAAAGGGTCTTACTGAAGTAACA
 CAAAGCTTAAAAATGGAAGATGAAGAGTTTAAAGAGAGTTTCAAGTATGCTACATCCAAAGCCCATCAGC
 TTGAGGAAGATATTGTGTCACTAACACATAAAGCAATTGAAAAAGAAACCGAATTAGACAGTTTAAAGGA
 CAAACTCAAGAAGGCACAACATGAAAGAGAACAACTTGAATGTCAGTTGAAGACAGAGAAGGATGAAAAG
 GAAGTTTATAAGGTACATTTGAAGAATACAGAAATAGAAAATACCAAGCTTATGTCAGAGGTCCAGACTT
 TAAAAAATTTAGATGGGAACAAAGAAAGCGTGATTACTCATTTCAAAGAAGAGATTGGCAGGCTGCAGTT
 ATGTTTGGCTGAAAAGGAAAAATCTGCAAGAAGCTTCTGCTTACAACCTCAAGTAAAGAAGATACTTGT
 TTTTAAAGGAGCAACTTCGTAAGCAGAGGAACAGGTTCAAGCAACTCGGCAAGAAGTTGTCTTTCTGG
 CTAAAGAACTCAGTGATGCTGTCAACGTACGAGACAGAACGATGGCAGACCTGCATACTGCACGCTTGA
 AAACGAGAAAGTGAAAAAGCAGTTAGCTGATGTCAGTGGCAGAACTTAACTAAATGCTATGAAAAAGAT
 CAGGACAAGACTGATACACTGGAACACGAACTAAGAAGAGAAGTTGAAGATCTGAAACTCCGTCTTCAGA
 TGGCTGCAGACCATTAATAAGAAAAATTAAGGAATGCCAAAGGCTCCAAAAACAAATAAACAACTTTC
 AGATCAATCAGCTAATAATAATAATGTCTTCACAAAGAAAACGGGGAATCAGCAGAAAGTGAATGATGCT
 TCAGTAAACACAGACCCAGCCACTTCTGCCTCTACTGTAGATGTAAGCCATCACCTTCTGCAGCAGAGG
 CAGATTTTGACATAGTAACAAAGGGCAAGTCTGTGAAATGACCAAGAAATGCTGACAAAACAGAAAA
 GTATAATAATGTAAACAACTCTTGAGGATGAGAAAGCAAAATGCAATAAATATGCTGATGAAGTTGCA
 AAAATGGAGCTGAAATGGAAGAACAAAGTGAAGTTGCTGAAATGTAAACTTGAAGTACTGAAGTAC
 AGGACAATTATAAAGAACTTAAAGGAGTCTAGAAAATCCAGCAGAAAGGAAATGGAAGATGGAGCAGA
 TGGTGCTTTTACCCAGATGAAATACAAAGGCCACCTGTCAGAGTCCCCTCTTGGGGACTGGAAGACAAT
 GTTGTCTGCAGCCAGCCTGCTCGAACTTTAGTCGGCCTGATGGCTTAGAGGACTCTGAGGATAGCAAG
 AAGATGAGAAATGTGCCTACTGCTCCTGATCCTCAAGTCAACATTTACGTGGGATGGGACAGGCTTTTG
 CTTTGATTCCAGCTTTGATGTTCAAGAAGTGTCCCCTCTGTGAGTTAATGTTTCTCCTAACTATGAT
 CAGAGCAAATTTGAAGACATGTTGAAAGTCACTGGAAGGTGTGCCCGATGTGCAGCGAGCAGTTCCCTC
 CTGACTATGACCAGCAGGTGTTTGAAGGCATGTGCAGACCCATTTTATCAGAATGTTCTAAATTTGA
 C

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234818 representing NM_001206901
Red=Cloning site Green=Tags(s)



MTS FQEVPLQTSNFAH VIFQNVAKSYLPNAHLECHYTLTPYIHPHPKDWVGIFKVGWSTARDYYTFLWSP
MPEHYVEGSTVNCVLA FQGYYLPNDDEGFYQFCYVTHKGEIRGASTPFQFRASSPVEELLTMEDEGNSDM
LVVTTKAGLLELKIEKTMKEKEELLKL IAVLEKETAQLREQVGRMERELNHEKERCDQLQAEQKGLTEVT
QSLKMENEEFKKRFSDATSKAHQLEEDIVSVTHKAI EKETELDSLKDKLKKAAQHEREQLEQCQLKTEKDEK
ELYKVHLKNT EIENTKLMSEVQTLKNLDGNKESVITHFKEEIGRLQLCLAEKENLQRTFLLTSSKEDTC
FLKEQLRKAAEEVQATRQEVVFLAKELSDAVNVDRDTMADLHTARLENEKVKKQLADAVAE LKLNAMKKD
QDKD TLEHEL RREVEDLKLRLQMAADHYKEKFKECQRLQKQINKLSDQSANNNNVFTKKTGNQKQVND A
SVNTDPATSASTVDVKPSPSAAEADFIVTKGQVCEMTKEIADKTEKYNKCKQLLQDEKAKCNKYADELA
KMELKQKEQVKIAENVGLELAEVDQNYKELKRSLPENARKMGMDGADGAFYPDEIQRPVVRVPSWGLEDN
VVCSSQARNF SRPDGLEDSEKDEDENVPATPPPSQHLRGMGTGFCFDDSSFVHKKCPLCELMFPFPNYD
QSKFEEHV SHWKVCPMCSEQFPDPDYDQOVFERHVOTHEDONVLFND

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 001206901

ORF Size: 2241 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001206901.1, NP_001193830.1

RefSeq Size: 3383 bp

RefSeq ORF: 2244 bp

Locus ID: 8887

UniProt ID: Q86VP1


Cytogenetics: 7p15.2

Protein Families: Druggable Genome

MW: 86.7 kDa

Gene Summary: This gene encodes a HTLV-1 tax1 binding protein. The encoded protein interacts with TNFAIP3, and inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degradation of this protein by caspase-3-like family proteins is associated with apoptosis induced by TNF. This protein may also have a role in the inhibition of inflammatory signaling pathways. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, May 2011]

Product images:



Circular map for RC234818