

Product datasheet for **RC234821**

KIFAP3 (NM_001204516) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIFAP3 (NM_001204516) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	KIFAP3
Synonyms:	dj190116.1; FLA3; KAP-1; KAP-3; KAP3; SMAP; Smg-GDS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234821 representing NM_001204516
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCAAGGGGAGGACGCCAGATACCTCAAAGCATTTCGACTTAAGAGTCTCAATGCCAACACAGATATAA
 CTTCCTGGCAAGGAAGGTGGTTGAAGAATGTAACTCATTTCCTTCAAACCTAAATGAGGTAGAACA
 GCTGTTGTACTATCTACAGAACCCTGATTCATTGTCAGGAAAAGAGAAAAAAGAAAAATCAAGCAAG
 CCTAAAGATCCACCTCCTTTGAAGGAATGGAGATTGATGAAGTTGCTAACATTAATGACATGGATGAAT
 ATATTGAGTTATTATGAAGATATTCCTGACAAAGTTCCGGGTTCTGCTTTGATCCTGCAGCTTCTCG
 AAATCCTGATAACTTGAAGAATACTATTGAATGAACTGCCCTTGGTGCATTAGCAAGGGTCTGAGA
 GAAGACTGGAAGCAAAGTGTGAGTTAGCTACAAACATAATTTACATCTTTTTTTGTTTCCAGCTTTT
 CTCATTTTCATGGACTTATTACTCACTATAAAATTGGAGCTCTGTGTATGAATATTATTGATCATGAGTT
 AAAAGACATGAGCTTTGGCAAGAAGAACTCTCAAAGAAGAAGAAAGCTGTTGATGAAGACCTGAAAAC
 CAAACCTTGAGAAAGGATTATGAAAAAACCTTTAAAAAGTACCAGGGGCTTGTGGTAAAAACGGAACAGC
 TATTACGAGTTGCTCTTTATTTGCTTCTGAATCTTGTCTGAGGATACTCGTACCGAACTGAAAATGAGGAA
 CAAGAACATAGTTCACATGTTGGTGAAGCCCTTGATCGGGACAATTTTGAGCTGCTAATTTTAGTTGTG
 TCATTTCTGAAGAACTCAGCATTTTTATGGAGAATAAAAATGATATGGTGGAAATGGATATTGTTGAAA
 AACTGGTGAAGATGATACCTTGTGAGCATGAAGACCTGCTGAATATCACCTCCGACTTTTACTAAACCT
 ATCCTTTGACACAGGACTGAGGAATAAGATGGTACAAGTTGGACTGCTTCCCAAGCTCACTGCACCTCTA
 GGCAATGACAACACAAACAAATAGCAATGTGTGTTCTTTACCACATAAGCATGGATGACCGCTTTAAAT
 CAATGTTTGCATACACTGACTGTATACCACAGTTAATGAAGATGCTGTTTGAATGTTGATGACACGAAT
 TGACTTGGAACTCAATTTCTTCTGCATTAATCTTGTCTGCTAACAAAAGAAATGTACAGCTTATCTGTGAA
 GGAAATGGGCTGAAGATGCTCATGAAGAGGCTCTGAAGTTTAAAGATCCATTGCTGATGAAAATGATTA
 GAAACATTTCTCAGCATGATGGACCAACTAAAAATCTGTTTATTGATTATGTTGGGGACCTTGCAGCCCA
 GATCTCTAATGATGAAGAAGAGGAGTTTGTGATTGAATGTTTGGGAACTTTGCAAACCTTGACCATTCCA
 GACTTAGACTGGGAATTGGTTCTTAAAGAATAAAGTTGGTTCCATACCTCAAGGATAAACTAAAACCGAG
 GTGCTGCAGAAGATGATCTTGTGTTTGAAGTGGTTATAATGATTGGAAGTGTATCCATGGATGACTCTTG
 TGCTGCATTGCTAGCCAAATCTGGCATAATCCCTGCACCTATTGAATTGCTAAATGCTCAACAAGAAGAT
 GATGAATTTGTGTGTCAGATAATTTATGTCTTCTACCAGATGGTTTTCCACCAAGCCACAAGAGACGTCA
 TAATCAAGGAAACACAGGCTCCAGCATATCTCATAGACCTAATGCATGATAAGAATAATGAAATCCGAAA
 GGTCTGTGATAATACATTAGATATTATAGCGGAATATGATGAAGAATGGGCTAAGAAAAATTCAGAGTGAA
 AAGTTTCGCTGGCATAACTCTCAGTGGCTGGAGATGGTAGAGAGTCGTCAGATGGATGAGAGTGAGCAGT
 ACTTGTATGGTGTGATCGAATTGAGCCATACATTCATGAAGGAGATATTCTCGAAAGACCTGACCTTTT
 CTACAACCTCAGATGGATTAATTGGCTCTGAAGGAGCCATAAGTCCCGATTCTTCAATGATTACCACCTT
 CAAAATGGAGATGTTGTTGGGCAGCATTCAATTTCTGGCAGCCTTGGAAATGGATGGCTTTGGCCAACCAG
 TTGGCATTCTTGGACGCCCTGCCACAGCATATGGATCCGCCCTGATGAACCTTACTACTATGGCTATGG
 ATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234821 representing NM_001204516
Red=Cloning site Green=Tags(s)


MQGEDARYLKSIRLKSLNANTDITSLARKVVEECKLIHPSKLNVEVQLLYLQNRRLDLSGKEKKEKSSK
 PKDPPPFEGMEIDEVANINDMDEYIELLYEDIPDKVRSALILQLARNPDNLEELLNETALGALARVLR
 EDWKQSVELATNIYIFFCFSSFSQFHGLITHYKIGALCMNIIDHELKRHELWQEELSFKKKAVDEDPEN
 QTLRKDYKTFKKYQGLVVKQEQLLRVALYLLNLAEDTRTELKMRNKNI VHMLVKALDRDNFELLILVV
 SFLKLLSIFMENKNDMVEMDIVEKLVKMIPCEHEDLLNITLRLLLNLSFDTGLRNKMVQVGLLPKLTALL
 GNDNYKQIAMCVLYHISMDDRFKSMFAYTDCIPQLMKMLFECSDERIDLELISFCINLAANKRNVLICE
 GNGLKMLMKRALKFKDPLLKMKIRNISQHDGPTKNLFDYVVDLAAQISNDEEEEFVIECLGTLANLTIP
 DLDWELVKEYKLVPLYKDKLKPAAEDDLVLEVVIMIGTVSMDDSCAALLAKSGIIPALIELLNAQQED
 DEFVCIYVYQVMVFHQATRDVIKETQAPAYLIDLMDKNNEIRKVCNTLDIIAEYDEEWAKKIQSE
 KFRWHNSQWLEMVESRQMDSEQYLYGDDRIEPIHEGDILERPDLFYNSDGLIASEGAI SPDFNFNDYHL
 QNGDVGQHSFPGSLGMDGFGQPVGILGRPATAYGFRPDEPYYYGYGS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001204516

ORF Size: 2244 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204516.1](#), [NP_001191445.1](#)

RefSeq Size: 2879 bp

RefSeq ORF: 2247 bp

Locus ID: 22920

UniProt ID: [Q92845](#)


Cytogenetics: 1q24.2

Protein Families: Druggable Genome

MW: 86.7 kDa

Gene Summary: The small G protein GDP dissociation stimulator (smg GDS) is a regulator protein having two activities on a group of small G proteins including the Rho and Rap1 family members and Ki-Ras; one is to stimulate their GDP/GTP exchange reactions, and the other is to inhibit their interactions with membranes. The protein encoded by this gene contains 9 'Armadillo' repeats and interacts with the smg GDS protein through these repeats. This protein, which is highly concentrated around the endoplasmic reticulum, is phosphorylated by v-src, and this phosphorylation reduces the affinity of the protein for smg GDS. It is thought that this protein serves as a linker between human chromosome-associated polypeptide (HCAP) and KIF3A/B, a kinesin superfamily protein in the nucleus, and that it plays a role in the interaction of chromosomes with an ATPase motor protein. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2011]

Product images:



Circular map for RC234821