

Product datasheet for **RC234837**

BRD2 (NM_001199456) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRD2 (NM_001199456) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BRD2
Synonyms:	BRD2-IT1; D6S113E; FSH; FSRG1; NAT; O27.1.1; RING3; RNF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC234837 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTTCGGTGCCTGCTTTGCAACTTACCCTGCCAACCCACCACCCCGGAGGTGTCCAATCCAAAA
 AGCCAGGACGAGTTACCAACCAGCTGCAATACCTACACAAGGTAGTGATGAAGGCTCTGTGGAAACATCA
 GTTCGCATGGCCATTCGGCAGCCTGTGGATGCTGCAAACTGGGTCTACCGGATTATCACAAAATTATA
 AAACAGCCTATGGACATGGGTACTATTAAGAGGAGACTTGAAAACAATTATTATTGGGCTGCTTCAGAGT
 GTATGCAAGATTTAATACCATGTTACCAACTGTTACATTTACAACAAGCCCACTGATGATATTGTCT
 AATGGCACAAACGCTGAAAAAGATATTCCTACAGAAGTTGCATCAATGCCACAAGAAGAACAAGAGCTG
 GTAGTGACCATCCCTAAGAACAGCCACAAGAAGGGGGCCAAGTTGGCAGCGCTCCAGGGCAGTGTACCA
 GTGCCCATCAGGTGCCTGCCGTCTTTCTGTGTACACACAGCCCTGTATACTCTCCACCTGAGATACC
 TACCACTGTCCTCAACATTCCCCACCCATCAGTCATTTCTCTCCACTTCTCAAGTCCTTGCCTGCTGCT
 GGACCCCGCTCCTTGCTGTTACTGCAGCTCCTCCAGCCAGCCCTTGCCAAGAAAAAGGCGTAAAGC
 GGAAAGCAGATACTACCACCCCTACACCTACAGCCATCTTGGCTCCTGGTCTCCAGCTAGCCCTCCTGG
 GAGTCTTGAGCCTAAGGCAGCACGGCTTCCCCCTATGCGTAGAGAGAGTGGTCGCCCCATCAAGCCCCCA
 CGCAAAGACTTGCTGACTCTCAGCAACAACACCAGAGCTCTAAGAAAGGAAAGCTTTCAGAACAGTTAA
 AACATTGCAATGGCATTTTGAAGGAGTACTCTCTAAGAAGCATGCTGCCTATGCTTGGCCTTTCTATAA
 ACCAGTGGATGCTTCTGCACTTGGCCTGCATGACTACCATGACATCATTAAAGACCCCATGGACCTCAGC
 ACTGTCAAGCGGAAGATGGAGAACCCTGATTACCGGGATGCACAGGAGTTTGTGCTGATGTACGGCTTA
 TGTCTCCAACCTGCTATAAGTACAATCCCCAGATCACGATGTTGTGGCAATGGCACGAAAGCTACAGGA
 TGTATTTGAGTTCGTTATGCCAAGATGCCAGATGAACCACTAGAACCAGGGCCTTTACCACTCTCTACT
 GCCATGCCCCCTGGCTTGGCCAAATCGTCTTACAGTCTCCAGTGAGGAAAGTAGCAGTGAAGCTCCT
 CTGAGGAAGAGGAGGAGGAAGATGAGGAGGACGAGGAGGAAGAAGAGAGTGAAAGCTCAGACTCAGAGGA
 AGAAAGGGCTCATCGCTTACGAGAACTACAGGAACAGCTTCGGGCAGTACATGAACAACCTGGCTGCTCTG
 TCCCAGGGTCCAATATCCAAGCCCAAGAGGAAAAGAGAGAAAAAGAGAAAAAGAAACGGAAGGCAG
 AGAAGCATCGAGGCCGAGCTGGGGCCGATGAAGATGACAAGGGGCTAGGGCACCCCGCCACCTCAACC
 TAAGAAGTCCAAGAAAGCAAGTGGCAGTGGGGTGGCAGTGTGCTTTAGGCCCTTCTGGCTTTGGACCT
 TCTGGAGGAAGTGGCACCAAGCTCCCCAAAAGGCCACAAGACAGCCCACTGCCCTGCCTACAGGTT
 ATGATTCAGAGGAGGAGGAAGAGAGCAGGCCATGAGTTACGATGAGAAGCGGCAGCTGAGCCTGGACAT
 CAACAAATTACCTGGGGAGAAGCTGGGCCGAGTTGTGCATATAATCCAAGCCAGGGAGCCCTCTTTACGT
 GATTCAAACCCAGAAGAGATTGAGATTGATTTTGAACACTCAAGCCATCCACACTTAGAGAGCTTGAGC
 GCTATGTCCTTTCTGCCTACGTAAGAAACCCCGAAGCCCTACACCATTAAGAAGCCTGTGGGAAAGAC
 AAAGGAGGAACCTGGCTTTGGAGAAAAAGCGGAATTAGAAAAGCGGTTACAAGATGTGAGCGGACAGCTC
 AATTCTACTAAAAAGCCCCCAAGAAAGCGAATGAGAAAACAGAGTCATCCTCTGCACAGCAAGTAGCAG
 TGTACGCCTTAGCGCTCCAGCTCCAGCTCAGATTCCAGCTCCTCTCTCTCTCGTCGTCGCTTTCAGA
 CACCAGTGATTCAGACTCAGGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234837 protein sequence
 Red=Cloning site Green=Tags(s)


MASVPALQLTPANPPPPEVSNPKKPGRVTNQLQYLHKVVMKALWKHQFAWPFRRQPVDAVKLGLPDYHKII
 KQPMDMGTIKRRELENNYYWAASECMQDFNTMFTNCYIYNKPTDDIVLMAQTLEKIFLQKVASMPQEEQEL
 VVTIPKNSHKKGAKLAALQGSVTSAHQVPVSSVSHALYTPPEIPTTVLNIPHPSVISSPLLKSLHSA
 GPPLLAVTAAPPAQPLAKKKGVKRAKADTTTPTPTAILAPGSPASPPGSLPEKAAARLPPMRRESGRPDKP
 RKDLPDSQQQHQSCKGKLEQLKHCNGILKELLSKKHAAYAWPFYKVPDASALGLHDYHDIKHPMDLS
 TVKRKMENTRDRDAQEFAADVRLMFSNCKYKYNPPDHDVVAMARKLQDVFEFRYAKMPDEPLEPGPLPVST
 AMPPLAKSSSSSSSESSSSSEEEEEDEDEDEDESESSDSEERAHRLAELQEQLRAVHEQLAAL
 SQGPISKPKRREKKEKKKRAEKHRGRAGADEDDKGPRAPPPQPKSKKASGSGGSAALGPSFGP
 SGGSGTKLPKKATKTAPPALPTGYDSESSSRPMSYDEKRLSLDINKLPGEKLRVVHIIQAREPSLR
 DSNPEEIEIDFETLKPSTLRELERVLSCLRKKPKPYTIKKPVGKTKELALEKKRELEKRLQDVSGQL
 NSTKPKPKANEKTESSSAQQVAVSRLSASSSSDSSSSSSSSSSSDTSDSDSG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001199456

ORF Size: 2262 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199456.1](#), [NP_001186385.1](#)

RefSeq Size: 3210 bp

RefSeq ORF: 2265 bp

Locus ID: 6046

UniProt ID: [P25440](#)


Cytogenetics: 6p21.32

Protein Families: Protein Kinase

MW: 83.2 kDa

Gene Summary: This gene encodes a transcriptional regulator that belongs to the BET (bromodomains and extra terminal domain) family of proteins. This protein associates with transcription complexes and with acetylated chromatin during mitosis, and it selectively binds to the acetylated lysine-12 residue of histone H4 via its two bromodomains. The gene maps to the major histocompatibility complex (MHC) class II region on chromosome 6p21.3, but sequence comparison suggests that the protein is not involved in the immune response. This gene has been implicated in juvenile myoclonic epilepsy, a common form of epilepsy that becomes apparent in adolescence. Multiple alternatively spliced variants have been described for this gene. [provided by RefSeq, Dec 2010]

Product images:



Circular map for RC234837