

Product datasheet for **RC234839**

WWP2 (NM_001270453) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WWP2 (NM_001270453) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	WWP2
Synonyms:	AIP2; WWp2-like
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234839 representing NM_001270453
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGCTGACCCTGAACCTGCAGACGGAGAACAAGGCAGCGTTGTCTCAGCGGAGAGCTGACAATTT
 TCCTGGACGGGCCAACTGTTGATCTGGGAAATGTGCCTAATGGCAGTGCCTGACAGATGGATCACAGCT
 GCCTTCGAGAGACTCCAGTGAACAGCAGTAGCTCCAGAGAACCAGGACCCAGCCCCAGCACAAATGC
 TTTGGTGAAGATCCCGGACGCACAGACATTGGGTGCTTCAGCCAGAACAACCCAGCAACCGGCGAGC
 AAAGCCCCGGTGTCTGGAGCCGGCACCAGCCCGTCAAGAAGTCCAGCCACAGTGGCTTGGCCAAATGG
 CACAGTGAATGATGAACCCACAACAGCCACTGATCCCGAAGAACCTTCCGTTGTTGGTGTGACGTCCCCA
 CCTGCTGCACCCTTGAGTGTGACCCCGAATCCCAACACGACTTCTCTCCCTGCCCCAGCCACACCGGCTG
 AAGGAGAGGAACCCAGCACTTCGGGTACACAGCAGCTCCAGCGGCTGCCAGGCCCCCGACGCTCTGCC
 TGCTGGATGGGAACAGCGAGAGCTGCCAACGGACGTGTCTATTATGTTGACCACAATACCAAGACCACC
 ACCTGGGAGCGGCCCTTCTCCAGGCTGGGAAAACGCACAGATCCCGAGGCAGGTTTTACTATGTGG
 ATCACAATACTCGGACCACCACCTGGCAGCGTCCGACCGCGGAGTACGTGCGCAACTATGAGCAGTGGCA
 GTCGCAGCGGAATCAGCTCCAGGGGGCCATGCAGCACTTCAGCCAAAGATTCTCTACCAGTCTTCGAGT
 GCTTCGACTGACCATGATCCCCTGGGCCCTCCCTCCTGGTGGGAGAAGAGACAGGACAATGGACGGG
 TGTATTACGTGAACCATAAACTCGCACGCCAGTGGGAGGATCCCGGACCCAGGGGATGATCCAGGA
 ACCAGCTCTGCCCCAGGATGGGAGATGAAATACACCAGCGAGGGGGTGCATACTTTGTGGACCACAAT
 ACCCGCACACCACCTTTAAGGATCCTCGCCCGGGTTTGAGTGGGGACGAAGCAAGTTCCCTTGGTG
 CTTATGACCGCAGTTTTCGGTGGAAGTATCACCAGTTCGTTTCTCTGCCATTCAAATGCCCTACCTAG
 CCACGTGAAGATCAGCGTTTCCAGGCAGACGCTTTTTCGAAGATTCTTCCAACAGATCATGAACATGAAA
 CCCTATGACCTGCGCCGCGGCTCTACATCATCATGCGTGGCGAGGAGGGCTGGACTATGGGGCATCG
 CCAGAGAGTGGTTTTTCTCCTGTCTCATGAGGTGCTCAACCCTATGTATTGTTTATTTGAATATGCCGG
 AAAGAACAATTACTGCCTGCAGATCAACCCCGCTCCTCCATCAACCCGACCCCTCACCTACTTTTCG
 TTTATAGGCAGATTCATCGCCATGGCGCTGTACCATGGAAAGTTCATCGACACGGGCTTCACCCTCCCTT
 TCTACAAGCGGATGCTCAATAAGAGACCAACCCTGAAAGACCTGGAGTCCATTGACCCTGAGTTCTACAA
 CTCCATTGTCTGGATCAAAGAGAACAACCTGGAAGAATGTGCCCTGGAGCTGTACTTCATCCAGGACATG
 GAGATACTGGGCAAGGTGACGACCCACGAGCTGAAGGAGGGCGGAGAGCATCCGGGTCACAGAGGAGA
 ACAAGGAAGAGTACATCATGCTGCTGACTGACTGGCGTTTCAACCGAGGGGTGGAAGAGCAGACCAAAGC
 CTTCTGGATGGCTTCAACGAGGTGGCCCCGCTGGAGTGGCTGCGCTACTTTGACGAGAAAGAGCTGGAG
 CTGATGCTGTGCGGCATGCAGGAGATAGACATGAGCGACTGGCAGAAGAGCACCATCTACCGGCACTACA
 CCAAGAACAGCAAGCAGATCCAGTGGTTCTGGCAGGTGGTGAAGGAGATGGACAACGAGAAGAGGATCCG
 GCTGCTGCAGTTTGTACCCGGTACCTGCCGCTGCCCGTGGGGGATTTGCCGAACCTCATCGGTAGCAAC
 GGACCACAGAAGTTTTGCATTGACAAAGTTGGCAAGGAAACCTGGCTGCCAGAAGCCACACCTGCTTCA
 ACCGTCTGGATCTCCACCCTACAAGAGCTACGAACAGCTGAGAGAGAAGCTGCTGTATGCCATTGAGGA
 GACCGAGGGCTTTGGACAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234839 representing NM_001270453
Red=Cloning site Green=Tags(s)


MQLTLNLQTENKGSVVSGGELTIFLDGPTVDLGNVPNGSALTDGSQLPSRDSSGTAVAPENRHQPPSTNC
 FGGRSRTHRHSGASARTTPATGEQSPGARSRHRQPVKNSGHSGLANGTVNDEPTTATDPEEPSVVGVTSP
 PAAPLSVTPNPNTTSLPAPATPAEGEEPSTSGTQQLPAAAQAPDALPAGWEQRELPNGRVYVVDHNTKTT
 TWERPLPPGWKERTDPRGRFYVDHNTRTTTWRPTAEYVRNYEQWQSQRNQLQGAMQHF SQRFLYQSSS
 ASTDHDPLGPLPPGWKRDNGRYYYYVNHNTRTTQWEDPRTQGMIQEPALPPGWEMKYTSEGVRYFVDHN
 TRTTTTFKDPRPGFESGKTQKQSPGAYDRSFRWKYHQRFLCHSNALPSHVKISVSRQTLFEDSFQQIMNMK
 PYDLRRRLYIIMRGEGLDYGGIAREWFFLLSHEVLNPMYCLFEYAGKNNYCLQINPASSINPDHLYFR
 FIGRFIAMALYHGKFI DTGFTLPFYKRLNKRPTLKDLESIDPEFYNSIVWIKENNL EECGLELYFIQDM
 EILGKVTTHLKEGGESIRVTEENKEEYIMLLTDWRFTRGVVEEQTKAFLDGFNEVAPLEWLRYFDEKELE
 LMLCGMQEIDMSDWQKSTIYRHYTKNSKIQWFVQVVKEMDNEKRIRLLQFVTGTCLRPVGGFAELIGSN
 GPQKFCIDKVGKETWLP RSHTCFNRLDLPPYKSYEQLREKLLYAIETEFGFGQE

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001270453

ORF Size: 2262 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001270453.1](#), [NP_001257382.1](#)

RefSeq Size: 4245 bp

RefSeq ORF: 2265 bp

Locus ID: 11060

UniProt ID: [O00308](#)

Cytogenetics: 16q22.1


Protein Families: Druggable Genome

Protein Pathways: Ubiquitin mediated proteolysis

MW: 86.8 kDa

Gene Summary: This gene encodes a member of the Nedd4 family of E3 ligases, which play an important role in protein ubiquitination. The encoded protein contains four WW domains and may play a role in multiple processes including chondrogenesis and the regulation of oncogenic signaling pathways via interactions with Smad proteins and the tumor suppressor PTEN. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 10. [provided by RefSeq, Jul 2012]

Product images:



Circular map for RC234839