

## Product datasheet for **RC234872**

### Metabotropic Glutamate Receptor 4 (GRM4) (NM\_001256813) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Metabotropic Glutamate Receptor 4 (GRM4) (NM_001256813) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Metabotropic Glutamate Receptor 4
Synonyms:	GPRC1D; mGlu4; MGLUR4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC234872 representing NM\_001256813  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCTGTGCCCTTGGGGCTCCTTGCTGGGCCTCTGCCCTGCCTGCCTGGGCTCCTCCCGCCTCCCC  
 ACAGATCACTGTTGACAAGGTTACTGAGTCAGCATGTAAAACCTGCAAAAATACCCAGATCAGCTACGC  
 CTCACAGCGCCAGACCTGAGTGACAACAGCCGCTACGACTTCTTCTCCCGCGTGGTGCCCTCGGACAG  
 TACCAGGCCAGGCCATGGTGGACATCGTCCGTGCCCTCAAGTGGAACTATGTGTCCACAGTGGCCTCGG  
 AGGGCAGCTATGGTGAGAGCGGTGTGGAGGCCCTCATCCAGAAGTCCCGTGGAGACGGGGCGTGTGCAT  
 CGCCAGTCGGTGAAGATACCACGGGAGCCCAAGGCAGGCAGTTCGACAAGATCATCCGCCCTCCTG  
 GAGACTTCGAACGCCAGGCAGTCATCATCTTTGCCAACGAGGATGACATCAGGCGTGTGCTGGAGGCAG  
 CACGAAGGGCCAACCAGACAGGCCATTTCTTCTGGATGGGCTCTGACAGCTGGGGCTCCAAGATTGCACC  
 TGTGCTGCACCTGGAGGAGGTGGCTGAGGGTGTGTACGATCCTCCCAAGAGGATGTCCGTACGAGGC  
 TTCGACCGCTACTTCTCCAGCCGACGCTGGACAACAACCGGCGCAACATCTGGTTTGCCGAGTTCTGGG  
 AGGACAATTCCACTGCAAGCTGAGCCGCCACGCCCTCAAGAAGGGCAGCCACGTCAAGAAGTGCACCAA  
 CCGTGAGCGAATTGGGCAGGATTCAGCTTATGAGCAGGAGGGGAAGGTGCAGTTTGTGATCGATGCCGTG  
 TACGCCATGGGCCACGCGCTGCACGCCATGCACCGTGACCTGTGTCCCGCCGCGTGGGGCTCTGCCCGC  
 GCATGGACCCTGTAGATGGCACCAGCTGCTTAAGTACATCCGAAACGTCAACTTCTCAGGCATCGCAGG  
 GAACCTGTGACCTTCAATGAGAATGGAGATGCGCCTGGGCGCTATGACATCTACCAATACCAGCTGCGC  
 AACGATTTGCCGAGTACAAGGTCATTGGCTCCTGGACTGACCACCTGCACCTTAGAATAGAGCGGATGC  
 ACTGGCCGGGAGCGGGCAGCAGCTGCCCGCTCCATCTGCAGCCTGCCAACCAGGGTGGAGCGGAA  
 GAAGACAGTGAAGGGCATGCCTTGCTGCTGGCACTGCGAGCCTTGACACAGGGTACCAGTACCAGGTGGAC  
 CGCTACACCTGTAAAGACGTGTCCCTATGACATGCGGCCACAGAGAACCACGCGGCTGCCGGCCATCC  
 CCATCATCAAGCTTGAGTGGGGCTCGCCCTGGGCGTGTGCCCTCTTCTGGCCGTGGTGGGCATCGC  
 TGCCACGTTGTTCTGTTGATCACCTTTGTGCGCTACAACGACACGCCCATCGTCAAGGCCCTCGGCCGT  
 GAACTGAGCTACGTGCTGCTGGCAGGCATCTTCTGTGCTATGCCACCACCTTCTCATGATCGCTGAGC  
 CCGACCTTGGCACCTGCTCGCTGCGCCGAATCTTCTGGGACTAGGGATGAGCATCAGCTATGCAGCCCT  
 GCTCACCAAGACCAACCGCATCTACCGCATCTTCGAGCAGGGCAAGCGCTCGGTCACTGCCCCACGCTTC  
 ATCAGCCCGCCTCACAGCTGGCCATCACCTTCAGCCTCATCTCGCTGCAGCTGCTGGGCATCTGTGTGT  
 GGTGTGTTGGTGGACCCCTCCCACTCGGTGGTGGACTTCCAGGACCAGCGGACACTCGACCCCGCTTCGC  
 CAGGGGTGTGCTCAAGTGTGACATCTCGGACCTGTGCTCATCTGCTGCTGGGTACAGCATGCTGCTC  
 ATGGTCACGTGCACCGTGTATGCCATCAAGACACGCGGCGTGCAGGACCTTCAATGAGGCCAAGCCCA  
 TTGGCTTACCATGTACACCACTTGATCGCTGGCTGGCCTTATCCCATCTTCTTTGGCACCTCGCA  
 GTCGGCCGACAAGCTGTACATCCAGACGACGACGCTGACGGTCTCGGTGAGTCTGAGCGCCTCGGTGTC  
 CTGGGAATGCTCTACATGCCAAAGTCTACATCATCTTCCACCCGGAGCAGAACGTGCCAAGCGCA  
 AGCGCAGCCTCAAAGCCGTGTTACGGCGGCCACCATGTCCAACAAGTTCACGCAGAAGGGCAACTCCG  
 GCCCAACGGAGAGGCCAAGTCTGAGCTCTGCGAGAACCTTGGGCCCCAGCGCTGGCCACCAACAGACT  
 TACGTCACCTACCAACCATGCAATC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC234872 representing NM\_001256813  
 Red=Cloning site Green=Tags(s)

```

MAVPLGAPC WASALPAWAPPGLPHRSLLTRLLSQHVKPAKIPQISYASTAPDLSDNSRYDFFSRVPSDT
YQAQAMVDIVRALKWNVSTVASEGSYGESGVEAFIQSREDGGVCI AQSVKIPREPKAGEFDKIIRLL
ETSNARAVIIFANEDDIRRVLEAARRANQTGHFFWMGSDSWGSKIAPVLHLEEVAEGAVTILPKRMSVRG
FDRYFSSRTL DNNRRNIWFAEFWEDNFHCKLSRHALKKGSHVKKCTNRERIGQDSAYEQEGKVQFVIDAV
YAMGHALHAMHRDLC PGRVGLCP RMDPVDGTQLLKYIRNVNFSGIAGNPVTFNENG DAPGRYDIYQYQLR
NDSAEYKVI GSWTDHLHLRIERMHWPGSGQQLPRSI CSLPCQPGERKKT VKGMPCCWHCEPCTGYQYQVD
RYTCKTCPYDMRPTENRTGCRPIPIIKLEWGSPWAVLPLFLAVVGIAATL FVVITFVRYNDTPIVKASGR
ELSYVLLAGIFLCYATTF LMI AEPDLGTC SLRRIFLGLGMSISYAALLTKTNRIYRIFEQ GKRSVSAPRF
ISPASQLAITFSLISLQLLGICVWFVVDPSHSVVD FQDQRTL DPRFARGVLKCDISDLSL ICLLGYSMLL
MVTCTVYAIKTRGVPET FNEAKPIGFTMYTTCIVWLA FIPIFFGTSQSADKLYIQTTTLTVSVLSASVS
LGMLYMPKVYIILFHPEQNVPKRKRSLKAVVTAATMSNKFTQKGNFRPNGEAKSEL CENLEAPALATKQT
YVTYTNHAI
  
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:**


**ACCN:** NM\_001256813

**ORF Size:** 2337 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001256813.3</a>
<b>RefSeq Size:</b>	6742 bp
<b>RefSeq ORF:</b>	2340 bp
<b>Locus ID:</b>	2914
<b>UniProt ID:</b>	<a href="#">Q14833</a>
<b>Cytogenetics:</b>	6p21.31
<b>Protein Families:</b>	Druggable Genome, GPCR, Transmembrane
<b>Protein Pathways:</b>	Neuroactive ligand-receptor interaction, Taste transduction
<b>MW:</b>	87.8 kDa
<b>Gene Summary:</b>	<p>L-glutamate is the major excitatory neurotransmitter in the central nervous system and activates both ionotropic and metabotropic glutamate receptors. Glutamatergic neurotransmission is involved in most aspects of normal brain function and can be perturbed in many neuropathologic conditions. The metabotropic glutamate receptors are a family of G protein-coupled receptors, that have been divided into 3 groups on the basis of sequence homology, putative signal transduction mechanisms, and pharmacologic properties. Group I includes GRM1 and GRM5 and these receptors have been shown to activate phospholipase C. Group II includes GRM2 and GRM3 while Group III includes GRM4, GRM6, GRM7 and GRM8. Group II and III receptors are linked to the inhibition of the cyclic AMP cascade but differ in their agonist selectivities. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]</p>

Product images:



Circular map for RC234872