

Product datasheet for **RC234954**

TAF1C (NM_001243156) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAF1C (NM_001243156) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TAF1C
Synonyms:	MGC:39976; SL1; TAFI95; TAFI110
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC234954 representing NM_001243156
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGACTTCCCCAGCTCCCTCCGCCTGCATTGTTTCTGACCGGCCCCCTTGGTCTGAGCGACGTCCTG
 ACCTCTCTTTCATGTGCAGCTGGCGAGACGCACTGACTCTGCCAGAGGCCAGCCCCAGAACTCAGAGAA
 TGGGGCACTGCATGTGACCAAGGACCTGCTGTGGGAGCCGGCAACCCCTGGGCTCTCCCATGCTGCCT
 CCCCTCATCGATCCCTGGGACCTGGCCTGACTGCCCGGGACCTGCTTTCCGCGGAGGGTACCGGTATC
 GGAAGCGGCCCGAGTCGTGCTGGATGTGACTGAGCAGATCAGCCGGTTCCTCTTGGATCATGGAGACGT
 AGCCTTTGCGCCCTGGGAAGCTGATGCTGGAGAATTTCAAGCTGGAGGGAGCGGGGAGCCGCACTAAG
 AAGAAGACAGTGGTCAGTGTGAAGAAGCTGCTCCAGGACCTCGGTGGACACCAGCCCTGGGGGTGTCCCT
 GGGCTTACCTCAGCAACCGACAGCGCCGTTCTCTATCCTCGGGGCCCATCTGGGCAGTCCGTGGC
 GAGCCACTTGGCAGAGCTGCTGCACGAGGAGCTGGTGTGCGGTGGGAGCAGCTGCTTCTGGATGAGGCC
 TGCACTGGGGCGCGCTGGCCTGGTTCCTGGAAGGACACCCAGTTCGGGCAGCTGGTCTACCCCTGCTG
 GAGGCGCCAGGACAGGCTGCATTTCAAGAGGTGCTTCTGACCCAGGTGACAATCCCCAATTCCTTGG
 GAAACCTGGACGCATCCAGCTCCAGGGACCTGTCCGGCAAGTGGTGCATGCACCGTCCAGGGAGAACT
 CTGCTGGCCGTCCTGACTACCACTGTGCCGTGTGGAAGTTTGGTAAACAGTGGCAGCCAACCCCTC
 TGCAGGCAATGCAGGTGGAGAAAGGGGCCACGGGGATCAGCCTCAGCCCTCACCTGCCGGGGAGCTGGC
 CATCTGCAGCCGCTCGGGAGCCGCTGCCTGTGGAGCCCTGAGGATGGGTGCGGCAAATCTACAGGGAC
 CCTGAGACCCTCGTGTCCGGGACTCCTCTCGTGGCGTTGGGCAGACTCACTGCGCACCCCTCGGGTGC
 TGACCGTGGGTGACCCGACCCGGAGTGAAGATGCTGGACACTCAGGGCCCGGGGCTGTGGTCTGTTGCT
 TTTTCGTTTGGGGCAGAGGCTTCGTGCCAGAAAGGGGAACGTGTCTGCTTACCCAGTACCTGGGGCAC
 TCCAGCCCCAAATGCCTCCCCCTACTCTTCATCTCGTCTGTACCCAGTCTCTCTCTACCTAGTGGACG
 AGCGCCTTCCCCTGGTCCGATGCTGAAGTGAACCATGGCTCCCCTCCCCGCTCTGCTGGCCCGACT
 GCTGCCTCCGCCCGGCCAGCTGCGTGCAGCCCTGCTCCTCGGAGGCCAGGTGGGCAGCTGCAGCTG
 CTGCACCTGGCAGGAGAAGGGCGTGGTGCCTGGCAGGCCCCCCCAGTCTCTTCTTCCAGGA
 TCGACTCCCTCCCTGCATTTCTCTGCTGGAGCCTAAGATCCAGTGGCGGCTGCAGGAGCGCCTGAAAGC
 ACCGACCATAGGTCTGGTGCCTGCTCCCGCCCTTGCCTCAGCGCCACACCAGGCTGGTGTCTTTC
 CAGCTCTCGGCGCGGGAGATGTCTTACCAGCAGCTCCGCCCCAGGTGGACTCCAGCCTCCGCAGAG
 ATGCTGGCCCTCCTGGCGACACCAACCTGACTGCCATGCCCCACAGCTTCTGGACCTCCAGGACAC
 TGCCGGCTGCAGCCAGTGGCTGAAGGCCCTGCTAAAAGTGCCCTGGCTCCTCCTGTGTGGACAGACCC
 ACCTTACCCACCGCCAGATGCTGGGCAGCACAGAGCTGCGGAGGGAGGAAGAGGAAGGGCAGCGCTGG
 GTGTGCTCCGCAAGGCCATGGCCCGAGGGCAGCTCCTGCTGCAGAGAGACCTGGGCTCCCTCCCTGCGGC
 AGAGCCACCCCTGCACCCGAGTCAGGCCTAGAGGACAAGCTCAGTGAAGCGCTGGGGGAAGCCTGGGCA
 GGCCGAGGGGCTGCCTGGTGGGAGAGGCAGGGCAGGACCTCGGAGCCCGGGAGACAGACCAGGGCGC
 CCAAGCGCCGACCCAGCTGTCCAGCAGCTTTTCGCTCAGTGGCCATGTGGATCCCTCAGAGGACACCAG
 CTCCCCATAGCCCTGAGTGGCCACCTGCTGATGCTTGCCTTGCCTTGCCTTGCCTTGCCTTGCCTTGCCT
 GAGTTGACTCCGGATGCATGCGCCAGGGCGTCCCATCAGAGCAGCGGCAGATGCTCCGTACTACATGG
 CCAAGCTACCACCCAGAGGGACACCCAGGCTGTGCCACACACCTCCCCTCCAGGCTCCAGCGT
 CCGGGCCACTCGCTCCAGCAGCACACCCGCTCTCTAGCTCTCAGCCCTCCGGAAGAAGCCTCGA
 ATGGGCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC234954 representing NM_001243156
 Red=Cloning site Green=Tags(s)

MDFPSSLRPALFLTGPLGLSDVDPDL S F M C S W R D A L T L P E A Q P Q N S E N G A L H V T K D L L W E P A T P G P L P M L P
 P L I D P W D P G L T A R D L L F R G G Y R Y R K R P R V V L D V T E Q I S R F L L D H G D V A F A P L G K L M L E N F K L E G A G S R T K
 K K T V V S V K K L L Q D L G G H Q P W G C P W A Y L S N R Q R R F S I L G G P I L G T S V A S H L A E L L H E E L V L R W E Q L L L D E A
 C T G G A L A W P G R T P Q F G Q L V Y P A G G A Q D R L H F Q E V V L T P G D N P Q F L G K P G R I Q L Q G P V R Q V V T C T V Q G E T
 L L A V R S D Y H C A V W K F G K Q W Q P T L L Q A M Q V E K G A T G I S L S P H L P G E L A I C S R S G A V C L W S P E D G L R Q I Y R D
 P E T L V F R D S S S W R W A D F T A H P R V L T V G D R T G V K M L D T Q G P P G C G L L L F R L G A E A S C Q K G E R V L L T Q Y L G H
 S S P K C L P P T L H L V C T Q F S L Y L V D E R L P L V P M L K W N H G L P S P L L L A R L L P P P R P S C V Q P L L L G G Q G G Q L Q L
 L H L A G E G A S V P R L A G P P Q S L P S R I D S L P A F P L L E P K I Q W R L Q E R L K A P T I G L A A V P P L S A P T P G L V L F
 Q L S A A G D V F Y Q Q L R P Q V D S S L R R D A G P P G D T Q P D C H A P T A S W T S Q D T A G C S Q W L K A L L K V P L A P P V W T A P
 T F T H R Q M L G S T E L R R E E E E G Q R L G V L R K A M A R G Q L L L Q R D L G S L P A A E P P P A P E S G L E D K L S E R L G E A W A
 G R G A A W W E R Q Q G R T S E P G R Q T R R P K R R T Q L S S S F S L S G H V D P S E D T S S P H S P E W P P A D A L P L P P T P P S Q
 E L T P D A C A Q G V P S E Q R Q M L R D Y M A K L P P Q R D T P G C A T T P P H S Q A S S V R A T R S Q Q H T P V L S S S Q P L R K K P R
 M G F

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:


Sgfl-MluI

Cloning Scheme:



ACCN:	NM_001243156
ORF Size:	2529 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243156.1 , NP_001230085.1
RefSeq Size:	3884 bp
RefSeq ORF:	2532 bp
Locus ID:	9013
UniProt ID:	Q15572
Cytogenetics:	16q24.1
Protein Families:	Transcription Factors
MW:	92.8 kDa
Gene Summary:	Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes the largest SL1-specific TAF. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2011]

Product images:



Circular map for RC234954