

Product datasheet for **RC234974**

Metabotropic Glutamate Receptor 4 (GRM4) (NM_001256811) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Metabotropic Glutamate Receptor 4 (GRM4) (NM_001256811) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Metabotropic Glutamate Receptor 4 |
| Synonyms: | GPRC1D; mGlu4; MGLUR4 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>RC234974 representing NM_001256811
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCTGGGAAGAGAGGCTTGGGCTGGTGGTGGGCCCGGTGCCCTTTGCCTGCTCCTCAGCCTTTACG
 GCCCTGGATGCCTTCTCCCTGGGAAAGCCAAAGGCCACCCTCACATGAATCCATCCGCATAGATGG
 GGACATCACACTGGGAGGCTGTTCCCGGTGCATGGCCGGGGCTCAGAGGGCAAGCCCTGTGGAGAACTT
 AAGAAGGAAAAGGGCATCCACCGGCTGGAGGCCATGCTGTTCCGCCTGGATCGCATCAACAACGACCCGG
 ACCTGCTGCCTAACATCACGCTGGGCGCCCGCATTCTGGACACCTGCTCCAGGGACACCCATGCCCTCGA
 GCAGTCGCTGACCTTTGTGACGGCGCTCATCGAGAAGGATGGCACAGAGGTCCGCTGTGGCAGTGGCGGC
 CCACCCATCATCACAAGCCTGAACGTGTGGTGGGTGTCATCGGTGCTTCAGGGAGCTCGGTCTCCATCA
 TGGTGGCCAAACATCCTTCGCCTTTCAAGATACCCAGATCAGCTACGCCTCCACAGCGCCAGACCTGAG
 TGACAACAGCCGCTACGACTTCTTCTCCGCGTGGTGCCTCGGACACGTACCAGGCCAGGCCATGGTG
 GACATCGTCCGTCGCCCTCAAGTGGAACTATGTGTCCACAGTGGCCTCGGAGGGCAGCTATGGTGAGAGCG
 GTGTGGAGGCCTTCATCCAGAAGTCCCCTGAGGACGGGGCGTGTGCATCGCCCAGTCCGGTGAAGATACC
 ACGGGAGCCAAAGGCAGGCGAGTTCGACAAGATCATCCGCCGCTCCTGGAGACTTCGAACGCCAGGGCA
 GTCATCATCTTTGCCAACGAGGATGACATCAGGCGTGTGCTGGAGGCAGCACGAAGGGCCAAACCAGACAG
 GCCATTTCTTCTGGATGGGCTCTGACAGCTGGGGCTCCAAGATTGCACCTGTGCTGCACCTGGAGGAGGT
 GGCTGAGGGTGTGTACGATCCTCCCCAAGAGGATGTCCGTACGAGACCGTGAGCGAATTGGGCAGGAT
 TCAGCTTATGAGCAGGAGGGGAAGGTGCAGTTTGTGATCGATGCCGTGTACGCCATGGGCCACGCGCTGC
 CCAGCTGCTTAAGTACATCCGAAACGTCAACTTCTCAGGATCGCAGGGAAACCCTGTGACCTTCAATGAG
 AATGGAGATGCGCCTGGGCGCTATGACATCTACCAATACCAGCTGCGCAACGATTCTGCCGAGTACAAGG
 TCATTGGCTCCTGGACTGACCACCTGCACCTTAGAATAGAGCGGATGCACTGGCCGGGAGCGGGCAGCA
 GCTGCCCGCTCCATCTGCAGCCTGCCCTGCCAACCGGTGAGCGGAAGAAGACAGTGAAGGGCATGCCCT
 TGCTGCTGGCACTGCGAGCCTTGACACAGGTTACCAGTACCAGGTGGACCGCTACACCTGTAAGACGTGTC
 CCTATGACATGCGGCCACAGAGAACCACAGGGCTGCCGGCCATCCCCATCATCAAGCTTGAGTGGGG
 CTCGCCCTGGGCGTGTGCCCTTCTCCTGGCCGTGGTGGGCATCGCTGCCACGTTGTTCTGTGGTGTGTC
 ACCTTTGTGCGCTACAACGACACGCCCATCGTCAAGGCCTCGGGCCGTGAAGTACGCTACGCTGCTGCTGG
 CAGGCATCTTCTGTGCTATGCCACCACCTTCTCATGATCGCTGAGCCCGACCTTGGCACCTGCTCGCT
 GCGCCGAATCTTCTGGGACTAGGGATGAGCATCAGTATGCAGCCCTGCTCACCAAGACCAACCGCATC
 TACCGCATCTTCGAGCAGGGCAAGCGCTCGGTGAGTGCACCGCTTTCATCAGCCCCGCTCACAGCTGG
 CCATCACCTTCAGCCTCATCTCGCTGCAGCTGTGGGCATCTGTGTGTGGTTTGTGGTGGACCCCTCCCA
 CTCGGTGGTGGACTTCCAGGACCAGCGGACACTCGACCCCGCTTCGCCAGGGGTGTGCTCAAGTGTGAC
 ATCTCGGACCTGTGCTCATCTGCCTGTGGGTACAGCATGTGCTCATGGTACGTGCACCGTGTATG
 CCATCAAGACACGCGGCTGCCGAGACCTTCAATGAGGCCAAGCCATTGGCTTACCATGTACACCAC
 TTGCATCGTCTGGCTGGCCTTCAATCCCATCTTCTTGGCACCTCGAGTCGGCCGACAAGCTGTACATC
 CAGACGACGACGCTGACGGTCTCGGTGAGTCTGAGCGCCTCGGTGCCCTGGGAATGCTCTACATGCCCA
 AAGTCTACATCATCTTCCACCCGGAGCAGAAGCTGCCCAAGCGCAAGCGCAGCCTCAAAGCCGTCGT
 TACGGCGCCACCATGTCCAACAAGTTCACGCAGAAGGGCAACTTCCGGCCCAACGGAGAGGGCAAGTCT
 GAGCTCTGCGAGAACCTTGGGCCCCAGCGCTGGCCACCAACAGACTTACGTCACTTACACCAACCATG
 CAATC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

| | |
|-------------------------------|--|
| ACCN: | NM_001256811 |
| ORF Size: | 2595 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001256811.2 , NP_001243740.1 |
| RefSeq Size: | 7293 bp |
| RefSeq ORF: | 2598 bp |
| Locus ID: | 2914 |
| UniProt ID: | Q14833 |
| Cytogenetics: | 6p21.31 |
| Protein Families: | Druggable Genome, GPCR, Transmembrane |
| Protein Pathways: | Neuroactive ligand-receptor interaction, Taste transduction |
| MW: | 96.6 kDa |
| Gene Summary: | L-glutamate is the major excitatory neurotransmitter in the central nervous system and activates both ionotropic and metabotropic glutamate receptors. Glutamatergic neurotransmission is involved in most aspects of normal brain function and can be perturbed in many neuropathologic conditions. The metabotropic glutamate receptors are a family of G protein-coupled receptors, that have been divided into 3 groups on the basis of sequence homology, putative signal transduction mechanisms, and pharmacologic properties. Group I includes GRM1 and GRM5 and these receptors have been shown to activate phospholipase C. Group II includes GRM2 and GRM3 while Group III includes GRM4, GRM6, GRM7 and GRM8. Group II and III receptors are linked to the inhibition of the cyclic AMP cascade but differ in their agonist selectivities. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012] |

Product images:



Circular map for RC234974