

Product datasheet for **RC235475A1V**

Human SNHG32 (NM_001287486) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNHG32 (NM_001287486) AAV Particle
Tag: Myc-DDK
Symbol: SNHG32
Synonyms: D6S57; G8
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC235475 representing NM_001287486
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGAGAGAAGCTTTGTATGGCTGTCATGCTTAGACAGTGATTCTGCAACTTGACCTTCAGGCTGGGAG
 AGGTGGAGAGCCATGCCTGTTCTCCTTCCTTGCATGGAATTTGCTGACACAATATCTCCGCCTGGTGC
 TGGGCATATCCTAAGAAGCTTACAACCTTCCTGTATTATCCTGTGTGAGCAGCTGTACCTTATTGGGGGA
 AAAATGCCTGAAAAT

ACGCGTACGCGCGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235475 representing NM_001287486
Red=Cloning site Green=Tags(s)

MERSFVWL SCLDSDSCNLTFR LGEVESHACSPSLLWNLLTQYLPPGAGHILRTYNFPVLSVSSCHLIGG
 KMPEN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001287486
ORF Size: 225 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001287486.1</u> , <u>NP_001274415.1</u>
RefSeq Size:	699 bp
RefSeq ORF:	228 bp
Locus ID:	50854
Cytogenetics:	6p21.33
MW:	8.3 kDa