

Product datasheet for **RC235477A1V**

Human SNHG32 (NM_001287488) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNHG32 (NM_001287488) AAV Particle
Tag: Myc-DDK
Symbol: SNHG32
Synonyms: D6S57; G8
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC235477 representing NM_001287488
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGAGAAGCTTTGTATGGCTGTCATGCTTAGACAGTGATTCTGCAACTTGACCTTCAGGCTGGGAG
AGGTGGAGAGCCATGCCTGTTCTCCTTCCTTGCTATGGAATTTGCTGACACAATATCTCCGCCTGGTGC
TGGGCATATCCTAAGAAGCTTACAACCTTCCTGTATTATCCTGTGTGAGCAGCTGCACCTATTGGGGGA
AAAATGCCTGAAAAT

ACGCGTACGCGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235477 representing NM_001287488
Red=Cloning site **Green**=Tags(s)

MERSFVWL SCLDSDSCNLTFR LGEVESHACSPSLLWNLLTQYLPPGAGHILRTYNFPVLSVSSCHLIGG
KMPEN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001287488
ORF Size: 225 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001287488.1</u> , <u>NP_001274417.1</u>
RefSeq Size:	496 bp
RefSeq ORF:	228 bp
Locus ID:	50854
Cytogenetics:	6p21.33
MW:	8.3 kDa