

Product datasheet for **RC235481**

PKIG (NM_001281445) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PKIG (NM_001281445) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PKIG
Synonyms: PKI-gamma
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235481 representing NM_001281445
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGATGGAGGTCGAGTCCTCTACTCGGACTTCATCTCCTGTGACCGGACAGGCCGTCGGAATGCGGTCC
CTGACATCCAGGGAGACTCAGAGGCTGTGAGCGTGAGGAAGCTGGCTGGAGACATGGCGAGCTGGCACT
CGAGGGGGCAGAAGGACAGGTGGAGGGAAGCGCCCGAGACAAGGAAGCTGGCAACCAGCCCGAGAGCAGC
GATGGGACCACCTCGTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235481 representing NM_001281445
Red=Cloning site Green=Tags(s)

MMEVESSYDFISCDRTGRRNAVPIQGDSEAVSVRKLADGMGELALEGAEGQVEGSAPDKEAGNQPQSS
DGTSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001281445

ORF Size: 228 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001281445.2
RefSeq Size:	1186 bp
RefSeq ORF:	231 bp
Locus ID:	11142
UniProt ID:	Q9Y2B9
Cytogenetics:	20q13.12
Protein Families:	Druggable Genome
MW:	8.4 kDa
Gene Summary:	This gene encodes a member of the protein kinase inhibitor family. Studies of a similar protein in mice suggest that this protein acts as a potent competitive cAMP-dependent protein kinase inhibitor, and is a predominant form of inhibitor in various tissues. The encoded protein may be involved in osteogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]