

Product datasheet for **RC235482**

SMIM8 (NM_001287445) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SMIM8 (NM_001287445) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SMIM8
Synonyms: C6orf162; dj102H19.2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235482 representing NM_001287445
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCTTCAGCACCTGAGCCTCCAACATTCAAAAAGGAACCCCAAAGAGAAAGAGTTTCAAAGCCCAG
GGCTCAGAGGGGTGCGCACAAACCTTATTTTCGTGCTGTGAATCCAGAGCTTTCATTAACCTAATGT
TTCTTATTCTCTGAAGCATCAGCATTTTGAAGTGTCTTTATGGAAGTGGGGACAGAAATGGGGTAAAA
ATGAGGTAAACAAAGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235482 representing NM_001287445
Red=Cloning site Green=Tags(s)

MSSAPEPPTFKKEPPKEKEFQSPGLRGVRTTTLFRAVNPELFIKPNVSYSLKHQHFEVVSFMEVGTEMGVK
MSVNKA


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001287445

ORF Size: 228 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|---|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001287445.2</u> |
| RefSeq Size: | 1306 bp |
| RefSeq ORF: | 231 bp |
| Locus ID: | 57150 |
| Cytogenetics: | 6q15 |
| Protein Families: | Transmembrane |
| MW: | 9 kDa |