

Product datasheet for **RC235520A1V**

Human Apc11 (ANAPC11) (NM_001289415) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Apc11 (ANAPC11) (NM_001289415) AAV Particle
Tag: Myc-DDK
Symbol: Apc11
Synonyms: APC11; Apc11p; HSPC214
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC235520 representing NM_001289415
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGGTGAAGATTAAGTGCTGGAACGGCGTGGCCACTTGGCTCTGGGTGGCCAACGATGAGAACTGTG
GCATCTGCAGGATGGCATTAAACGGATGCTGCCCTGACTGCAAGGTGCCCGGCGACGACTGCCCGCTGGT
GTGGGGCCAGTGCTCCCACTGCTCCACATGCATTGCATCTCAAGTGGCTGCACGCACAGCAGGTGCAG
CAGCACTGCCCATGTGCCGCCAGGAATGGAAGTCAAGGAG

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235520 representing NM_001289415
Red=Cloning site Green=Tags(s)

MKVKIKCWNGVATWLWVANDENCGICRMAFNCCPDCKVPGDDCPLVWGQCSHCFMHCIKWLHAQQVQ
QHCPMCRQEWKFKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001289415
ORF Size: 252 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_001289415.1](#), [NP_001276344.1](#)

RefSeq Size: 885 bp

RefSeq ORF: 255 bp

Locus ID: 51529

UniProt ID: [Q9NYG5](#)

Cytogenetics: 17q25.3

MW: 9.8 kDa