

## Product datasheet for **RC235531**

### TSC22D1 (NM\_001243798) Human Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
Product Name: TSC22D1 (NM\_001243798) Human Tagged ORF Clone

Tag: Myc-DDK

Symbol: TSC22D1

Synonyms: Ptg-2; TGFB114; TSC22

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Cell Selection: Neomycin

ORF Nucleotide Sequence: >RC235531 representing NM\_001243798  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGATCTAGTAAAAGCCATTTGATGTATGCGGTCAGAGAAGAAGTGGAGGTCCTCAAAGAGCAAATCA  
AAGAATAATAGAGAAAATTCCAGCTGGAGCAGGAGAACAACTGCTGAAGACTGGCCAGTCTCTGA  
GCAGCTTGCCAGTTTCAGGCCAGCTGCAGACTGGCTCCCCCTGCCACCACCCAGCCACAGGGCACC  
ACACAGCCCCCGCCAGCCAGCATCGCAGGGCTCAGGACCAACCGCA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235531 representing NM\_001243798  
Red=Cloning site Green=Tags(s)

MDLVKSHLMYAVREEVEVLKEQIKELIEKNSQLEQENLLKTLASPEQLAQFQAQLQTGSPATTQPQGT  
TQPPAQPASQSGPTA

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001243798

ORF Size: 258 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001243798.1, NP_001230727.1</u>
<b>RefSeq Size:</b>	2877 bp
<b>RefSeq ORF:</b>	261 bp
<b>Locus ID:</b>	8848
<b>Cytogenetics:</b>	13q14.11
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	9.8 kDa
<b>Gene Summary:</b>	This gene encodes a member of the TSC22 domain family of leucine zipper transcription factors. The encoded protein is stimulated by transforming growth factor beta, and regulates the transcription of multiple genes including C-type natriuretic peptide. The encoded protein may play a critical role in tumor suppression through the induction of cancer cell apoptosis, and a single nucleotide polymorphism in the promoter of this gene has been associated with diabetic nephropathy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2011]