

Product datasheet for **RC235622**

CCDC103 (NM_001258398) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CCDC103 (NM_001258398) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: CCDC103
Synonyms: CILD17; PR46b; SMH
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235622 representing NM_001258398
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAAGAATGACATCATCACTTCAAGGCTTTGGAGAAAGAGCTGCAGGCTGCACTCACTGCTGATG
AGAAGTACAAACGGGAGAATGCTGCCAAGTTACGGGCAGTGGAACAGAGGGTGGCTTCTATGAGGAGTT
CAGGGGTATTGTCCTTGCATCACATCTGAAGCCACTGGAGCGGAAGGATAAGATGGGAGGAAAGAGAACT
GTGCCCTGGAAGTGTCACTATTTCAGGAAGGACCTTCCAGGATGTGGCCACTGAAATCTCCCCGGTAG
AAAGCCCCCTCCAGCCGAGACGTCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235622 representing NM_001258398
Red=Cloning site Green=Tags(s)

MERNDIINFKALEKELQAALTADEKYKRENAAKLRAVEQRVASYEEFRGIVLASHLKPLERKDKMGGKRT
VPWNCHTIQGRTFQDVATEISPVESPPARDVC

TRTRPLEQK**L**ISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001258398

ORF Size: 309 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001258398.2</u>
RefSeq Size:	1751 bp
RefSeq ORF:	312 bp
Locus ID:	388389
Cytogenetics:	17q21.31
MW:	12.2 kDa
Gene Summary:	This gene encodes a protein that contains a coiled-coil domain. [provided by RefSeq, Apr 2012]