

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for RC235643

## Cyclophilin A (PPIA) (NM\_001300981) Human Tagged ORF Clone

# **Product data:**

Product Type:	Expression Plasmids
Product Name:	Cyclophilin A (PPIA) (NM_001300981) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PPIA
Synonyms:	CYPA; CYPH; HEL-S-69p
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	<pre>&gt;RC235643 representing NM_001300981 Red=Cloning site Blue=ORF Green=Tags(s)</pre>
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C
	ATGTGTCAGGGTGGTGACTTCACACGCCATAATGGCACTGGTGGCAAGTCCATCTATGGGGAGAAATTTG AAGATGAGAACTTCATCCTAAAGCATACGGGTCCTGGCATCTTGTCCATGGCAAATGCTGGACCCAACAC AAATGGTTCCCAGTTTTTCATCTGCACTGCCAAGACTGAGTGGGTTGGATGGCAAGCATGTGGTGTTTGGC AAAGTGAAAGAAGGCATGAATATTGTGGAGGCCATGGAGCGCTTTGGGTCCAGGAATGGCAAGACCAGCA AGAAGATCACCATTGCTGACTGTGGACAACTCGAA
	ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAG <b>GTTTAA</b>
Protein Sequence:	>RC235643 representing NM_001300981 <mark>Red</mark> =Cloning site Green=Tags(s)
	MCQGGDFTRHNGTGGKSIYGEKFEDENFILKHTGPGILSMANAGPNTNGSQFFICTAKTEWLDGKHVVFG KVKEGMNIVEAMERFGSRNGKTSKKITIADCGQLE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	Sgfl-Mlul



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US




#### **Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

## Plasmid Map:



ACCN:	
ORF Size:	
OTI Disclaimer:	

## NM\_001300981

## 315 bp

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

<b>Cyclophilin A (PPIA) (NM_001300981) Human Tagged ORF Clone – RC235643</b>	
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol> <li>Centrifuge at 5,000xg for 5min.</li> <li>Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>Close the tube and incubate for 10 minutes at room temperature.</li> <li>Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
RefSeq:	<u>NM 001300981.2</u>
RefSeq Size:	2445 bp
RefSeq ORF:	318 bp
Locus ID:	5478
UniProt ID:	<u>P62937</u>
Cytogenetics:	7p13
MW:	11.9 kDa
Gene Summary:	This gene encodes a member of the peptidyl-prolyl cis-trans isomerase (PPlase) family. PPlases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. The encoded protein is a cyclosporin binding-protein and may play a role in cyclosporin A-mediated immunosuppression. The protein can also interact with several HIV proteins, including p55 gag, Vpr, and capsid protein, and has been

shown to be necessary for the formation of infectious HIV virions. Multiple pseudogenes that

map to different chromosomes have been reported. [provided by RefSeq, Jul 2008]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US