

Product datasheet for **RC235677**

MAGP2 (MFAP5) (NM_001297712) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MAGP2 (MFAP5) (NM_001297712) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: MFAP5
Synonyms: AAT9; MAGP-2; MAGP2; MFAP-5; MP25
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235677 representing NM_001297712
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCGCTCTTGGACCCAAGGTGCTGCTGTTTCTTGCTGCATTCATCATCACCTCTGACTGGATACCC
TGGGGTCAATAGTCAACGAGGAGACGATGTGACTCAAGCGACTCCAGAAACATTCACAGAAGATCCTAA
TCTGGTGAATGATCCCGCTACAGATGAAACAGTTTGGCTGTTTGGCTGATATTGCACCTCCACAGAT
GACTTGGATGAGCTTTGCCGTCAGATGGCTGGCTGCCCCCTAGGAGACTCCGTCGCTCCAATTACTTCC
GACTTCTCCCTGTGAAAATGTGGATTGCAGAGACCCAATGGTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235677 representing NM_001297712
Red=Cloning site Green=Tags(s)

MSLLGPKVLLFLAAFIIITSDWIPLGVNSQRGDDVTQATPETFTEDPNLVNDPATDETVLAVLADIAPSTD
DLDELCRQMAGLPPRRLRRSNYFRLPPCENVDLQRPNGL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001297712

ORF Size: 327 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001297712.2</u>
RefSeq Size:	2757 bp
RefSeq ORF:	330 bp
Locus ID:	8076
Cytogenetics:	12p13.31
Protein Families:	Secreted Protein
MW:	12.5 kDa
Gene Summary:	This gene encodes a 25-kD microfibril-associated glycoprotein which is a component of microfibrils of the extracellular matrix. The encoded protein promotes attachment of cells to microfibrils via alpha-V-beta-3 integrin. Deficiency of this gene in mice results in neutropenia. Alternate splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014]