

Product datasheet for RC235701

COMMD4 (NM_001284378) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: COMMD4 (NM_001284378) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: COMMD4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC235701 representing NM_001284378
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGGTTCCGGTCTGTGGTGATCTGGACTGTCCCGACTGGGTCCTGGCAGAAATCAGCACGCTGGCCA
 AGATGTCCTCTGTGAAGTTGCGGCTGCTCTGCAGCCAGGTAATAAGGAGCTGCTGGGACAGGGGATTGA
 TTATGAGAAGATCCTGAAGCTCACGGCTGACGCCAAGTTTGAGTCAGGCGATGTGAAGGCCACAGTGGCA
 GTGCTGAGTTTCATCTCTCCAGTGC GGCCAAGCACAGTGTGATGGCGAATCCTTGTCAGTGAAGTGC
 AGCAGCTGGGGCTGCCAAAGAACTGAAGCAGGCCAGACCCTGATGAGCTCCCTGGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC235701 representing NM_001284378
 Red=Cloning site Green=Tags(s)

MRFRFCGDLDCPDWVLAELSTLAKMSSVKLRLLCSQVLKELLGQGIDYEKILKLTADAKFESGDVKATVA
 VLSFILSSAAKHSVDGESLSSELQQLGLPKELKQAQTLMSLGL

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_001284378

ORF Size: 339 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001284378.2](#)

RefSeq Size:	708 bp
RefSeq ORF:	342 bp
Locus ID:	54939
UniProt ID:	Q9H0A8
Cytogenetics:	15q24.2
MW:	12.7 kDa
Gene Summary:	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). Down-regulates activation of NF-kappa-B.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC235701