

## Product datasheet for **RC235712**

### ZCCHC10 (NM\_001300822) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZCCHC10 (NM_001300822) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZCCHC10
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC235712 representing NM_001300822 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCGACTCCCATGCATCGGCTAATAGCCCGGAGACAAGCATTGCACACAGAGTTGCAGCCTGTCAAGA  
CCTTTTGGATCCTGATTCAGCCATCCATCGTTATTAGTGAAGCAAATAAGCAACATGTAAGATGTCAGAA  
ATGCTTGGAAATTTGGACATTGGACTTATGAATGCACAGGAAAAAGAAAAATACCTACATAGGCCCTCAAGG  
ACAGCAGAACTAAAGAAAGCTTTAAAAGAAAAAGAAAACAGATTATTATTGCAACAAAGGAGCTTTTTTC  
CCCTCGTGTTCATCAGCATTGGAGAAACCAATGTAGAAAGAAAGGCCAAGAAAAAAGGTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC235712 representing NM_001300822 Red=Cloning site Green=Tags(s)
-------------------	---

MATPMHRLIARRQAFDELQPVKTFWILIQPSIVISEANKQHVRQCQKCLEFGHWTYECTGKRKYLHRPSR  
TAEKKALKKEKENRLLLQQRSFPPRVYQHWRNQCRRKKGQEKV

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	Sgfl-MluI
--------------------	-----------




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001300822  
 ORF Size: 342 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001300822.3</a>
<b>RefSeq Size:</b>	2286 bp
<b>RefSeq ORF:</b>	345 bp
<b>Locus ID:</b>	54819
<b>UniProt ID:</b>	<a href="#">Q8TBK6</a>
<b>Cytogenetics:</b>	5q31.1
<b>MW:</b>	14.3 kDa