

Product datasheet for **RC235721**

C1orf150 (GCSAML) (NM_001281835) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C1orf150 (GCSAML) (NM_001281835) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: GCSAML
Synonyms: C1orf150
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235721 representing NM_001281835
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

**ATGGGAAATTATCTCCTGCGAAAACCTCAGGCAGGAAATGACTACATTTGAAAGAAAACCTCAAGATCAAG
ATAAGAAAAGCCAAGAAGTTTCATCCACTTCTAATCAGGAAAACGAGAATGGCAGTGGTTCTGAAGAAGT
GTGCTACACTGTCATTAATCACATCCCCATCAGAGATCCTCCCTGAGCTCCAATGATGATGGCTATGAG
AACATTGACTCCCTCACAAAGGAAAGTGAGACAGTTTAGAGAAAGGTCAGAGACAGAATATGCCCTTCTTA
GGACTTCTGTTAGTAGGCCTTGTCTGCACCCATGAGCATGATTATGAAGTTGTGTTTCCACAC**

**ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**

Protein Sequence: >RC235721 representing NM_001281835
Red=Cloning site Green=Tags(s)

MGNLLRKLRLQEMTTFERKLQDQDKKSQEVSSSTSNQENENSGSSEEVICYTVINHIPHQRSSLSSNDGGE
NIDSLTRKVRQFRERSETEYALLRTSVSRPCSCTHEHDYEVVFPH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001281835

ORF Size: 345 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001281835.2</u>
RefSeq Size:	4094 bp
RefSeq ORF:	348 bp
Locus ID:	148823
UniProt ID:	<u>Q5JQS6</u>
Cytogenetics:	1q44
MW:	13.9 kDa
Gene Summary:	This gene encodes a protein thought to be a signaling molecule associated with germinal centers, the sites of proliferation and differentiation of mature B lymphocytes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]