

Product datasheet for **RC235784A1V**

Human C11orf51 (ANAPC15) (NM_001278491) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C11orf51 (ANAPC15) (NM_001278491) AAV Particle
Tag: Myc-DDK
Symbol: C11orf51
Synonyms: APC15; C11orf51; HSPC020
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC235784 representing NM_001278491
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGTCCACTTTGTTCCCCTCACTCTTCCTCGTGTGACTGAGACTCTGTGGTTAATCTGGATCGACCCT
 GTGTGGAAGAGACAGAGCTGCAGCAGCAGGAACAGCAGCATCAGGCCTGGCTCCAAGCATCGCGGAGAA
 AGACAACAACCTGGTTCTATTGGCAAGCCAGCCTCAGAGCACTATGATGACGAGGAAGAAGAGGATGAT
 GAAGATGATGAGGATAGTGAAGAGGACTCAGAGGATGATGAGGATATGCAGGACATGGACGAGATGAATG
 ACTACAATGAGTCACCGGATGATGGAGAGGTCAATGAGGTGGACATGGAAGGCAACGAACAGGATCAGGA
 CCAGTGGATGATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235784 representing NM_001278491
Red=Cloning site Green=Tags(s)

MSTLFPSLFPRVTETLWFNLDRPCVEETELQQQEQQHQAWLQSI AEKDNLLVPIGKPASEHYDDEEEEDD
 EDDEDESEDEDEDMDDEMNDYNESPDDGEVNEVDMEGNEQDQDQWMI

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001278491



[View online »](#)

ORF Size:	363 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001278491.1</u> , <u>NP_001265420.1</u>
RefSeq Size:	780 bp
RefSeq ORF:	366 bp
Locus ID:	25906
UniProt ID:	<u>P60006</u>
Cytogenetics:	11q13.4
MW:	14.3 kDa