

Product datasheet for **RC235808**

Sigma1 receptor (SIGMAR1) (NM_001282206) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sigma1 receptor (SIGMAR1) (NM_001282206) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SIGMAR1
Synonyms: ALS16; DSMA2; hSigmaR1; OPRS1; SIG-1R; sigma1R; SR-BP; SR-BP1; SRBP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235808 representing NM_001282206
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCGGGTGGCTGGATGGGCGCCATGTGCCTTCTGCACGCCTCGCTGTCCGAGGCGCTACTGGGCTGAGA
TCTCGGATACCATCATCTCTGGCACCTCCACCAGTGGAGAGAGGGCACCACAAAAGTGAGGTCTTCTA
CCCAGGGGAGACGGTAGTACACGGGCCTGGTGAGGCAACAGCTGTGGAGTGGGGCCAAACATGGATG
GTGGAGTACGGCCGGGGGTCATCCCATCCACCCTGGCCTTCGCCTGGCCGACACTGTCTTCAGCACCC
AGGACTTCTCACCTTCTATACTCTTCGCTCTATGCTCGGGCCTCCGGCTTGAGCTCACACCTA
CCTCTTTGGCCAGGACCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235808 representing NM_001282206
 Red=Cloning site Green=Tags(s)

MRVAGWAPCACTPRCPRRYWAEISDTIISGTFHQWREGTTKSEVFYPGETVVHGPGEATAVEWGPNTWM
VEYGRGVIPSTLAFALADTVFSTQDFLTLFYTLRSYARGLRLELTTYLFGQDP

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001282206

ORF Size: 369 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282206.1 , NP_001269135.1
RefSeq Size:	1681 bp
RefSeq ORF:	372 bp
Locus ID:	10280
UniProt ID:	Q99720
Cytogenetics:	9p13.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	14.3 kDa
Gene Summary:	This gene encodes a receptor protein that interacts with a variety of psychotomimetic drugs, including cocaine and amphetamines. The receptor is believed to play an important role in the cellular functions of various tissues associated with the endocrine, immune, and nervous systems. As indicated by its previous name, opioid receptor sigma 1 (OPRS1), the product of this gene was erroneously thought to function as an opioid receptor; it is now thought to be a non-opioid receptor. Mutations in this gene has been associated with juvenile amyotrophic lateral sclerosis 16. Alternative splicing of this gene results in transcript variants encoding distinct isoforms. [provided by RefSeq, Aug 2013]