

Product datasheet for **RC235820**

C2orf43 (LDAH) (NM_001282724) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C2orf43 (LDAH) (NM_001282724) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C2orf43
Synonyms: C2orf43; hLDAH
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC235820 representing NM_001282724
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGACTCAGAACTCAAGGAAGAAATTCCTGTGCATGAGGAATTCATTTTGTGTGGTGGAGCCGAAACCC
AGGTTCTAAAATGTGGGCCCTGGACAGACCTCTTTCATGATCAAAGTGCAAAAGGCCTAAGCTGCTTAT
TTTCATTATTCCTGGTGAGCGTAAAGTCTGGATTCAAACGCTCAAGAAATTAAGGACATTTATGGACTA
AATGGACAAATAGAGCACAACTAGCTTTCCTGAGAACTCATGTGCCAAAGGACATGAACTTGTGCTCA
TTGGCCATTCAATAGGCAGCTATTTCACTTCAGATGCTGAAGCGAGTCCCTGAGCTCCCGCTAATGCT
GCCTACCTTGGGGCCAAGAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235820 representing NM_001282724
Red=Cloning site Green=Tags(s)

MDSELKEEIPVHEEFILCGGAETQVLKCGPWTDLFHDQSVKRPKLLIFIIPGERKSLDSNAQEIKDIYGL
 NGQIEHKLAFLRTHVPKDMKLVLIHSGISYFTLQMLKRVPELPLMLPTLGAKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_001282724

ORF Size: 372 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001282724.2](#)

RefSeq Size: 3602 bp

RefSeq ORF: 375 bp
 Locus ID: 60526
 Cytogenetics: 2p24.1
 Protein Families: Druggable Genome
 MW: 14.5 kDa

Gene Summary: Serine lipid hydrolase associated with lipid droplets. Highly expressed in macrophage-rich areas in atherosclerotic lesions, suggesting that it could promote cholesterol ester turnover in macrophages.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC235820