

Product datasheet for **RC235826A1V**

Human LINC02901 (NM_001302839) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LINC02901 (NM_001302839) AAV Particle
Tag: Myc-DDK
Symbol: LINC02901
Synonyms: yR211F11.1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC235826 representing NM_001302839
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGAACGAGGAACACTGCCCGGGGAGTTTCTTCTCTGCAAATAAGGGAGTGTGTACTGAATTACC
 GTTCCAGCTCCAACATCCAGGATTTCCAGCACTACCTGCAGAGTAGTGGGAGGAGAGATAGAGGAAGAAG
 TGAAGATAAGAAACCACTCGAAGCCGGTGTCTGGTGTGGGACCGTGGTGGTGGGACCGCAGCAGTCGT
 GCTGTCCATCTGCTCTCCGAGGGGTGGCACATCCTAGCCTCTATCTTTCCAGAGAAGATCCCCCTC
 GTCTGCTTTCCCGTCTCAGCCTTTTAGTGTGTGAACAGTTCTGGTGTACTCAGCAACGCTTTTGTT
 AGCGCCTTTACCAGCAAGTACTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAAGTTTAA

Protein Sequence: >RC235826 representing NM_001302839
 Red=Cloning site Green=Tags(s)

MKNEEHCPGSFFLCKIRECVLNYRFQLQHPGFQHYLQSSGRRDRGRSEDKKPLEAGVWCWDRGGWDGSSR
 AVHLLFRGVAHPSLYLFPREDPPRLLFPRLSLLVCEQFWCYSATLLLAPLPASTC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001302839



[View online »](#)

ORF Size:	375 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001302839.1</u> , <u>NP_001289768.1</u>
RefSeq Size:	988 bp
RefSeq ORF:	378 bp
Locus ID:	100130967
UniProt ID:	<u>Q4VX62</u>
Cytogenetics:	6q25.3
MW:	14.5 kDa