

## Product datasheet for **RC235831**

### **C11orf51 (ANAPC15) (NM\_001278486) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C11orf51 (ANAPC15) (NM\_001278486) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ANAPC15  
**Synonyms:** APC15; C11orf51; HSPC020  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC235831 representing NM\_001278486  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

**ATGTCCACTTTGTTCCCTCACTCTTCCCTCGTGTGACTGAGACTCTGTGGTTTAACTGGATCGACCT**  
**GTGTGGAAGAGACAGAGCTGCAGCAGCAGGAACAGCAGCATCAGGCTGGCTCCAAGCATCGCGGAGAA**  
**AGACAACAACCTGGTTCCTATTGCAAGCCAGCCTCAGAGCACTATGATGACGAGGAAGAAGAGGATGAT**  
**GAAGATGATGAGGATAGTGAAGAGGACTCAGAGGATGATGAGGATATGCAGGACATGGACGAGATGAATG**  
**ACTACAATGAGTCACCGGATGATGGAGAGGTCAATGAGCCTCCTCTCCAGGTGGACATGGAAGGCAACGA**  
**ACAGGATCAGGACCAGTGGATGATC**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC235831 representing NM\_001278486  
**Red**=Cloning site **Green**=Tags(s)

MSTLFPSLFPRVTETLWFNLDRPCVEETELQQQEQQHQAWLQSI AEKDNLLVPIGKPASEHYDDEEEEDD  
EDDEDSEEDSEDDMDEMNDYNESPDDGEVNEPPLQVDMEGNEQDQDQWMI


**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


ACCN: NM\_001278486

ORF Size: 375 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001278486.2</a>
<b>RefSeq Size:</b>	889 bp
<b>RefSeq ORF:</b>	378 bp
<b>Locus ID:</b>	25906
<b>UniProt ID:</b>	<a href="#">P60006</a>
<b>Cytogenetics:</b>	11q13.4
<b>MW:</b>	15.2 kDa
<b>Gene Summary:</b>	Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle. In the complex, plays a role in the release of the mitotic checkpoint complex (MCC) from the APC/C: not required for APC/C activity itself, but promotes the turnover of CDC20 and MCC on the APC/C, thereby participating in the responsiveness of the spindle assembly checkpoint. Also required for degradation of CDC20.[UniProtKB/Swiss-Prot Function]