

Product datasheet for **RC235838A1V**

Human POLR1H (NM_001278785) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human POLR1H (NM_001278785) AAV Particle
Tag:	Myc-DDK
Symbol:	POLR1H
Synonyms:	A12.2; HTEX-6; HTEX6; hZR14; Rpa12; tctex-6; TCTEX6; TEX6; ZNRD1; ZR14
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC235838 representing NM_001278785 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC ATGCTGTGCATGGACCTCGCCAATACTTGCTCCAGCTTTCAGTCGGACCTGGATTTCTGTTCAGATTGCGCTCGGTCCTGCCTCTGCCGGGGCTCAGGATACGGTCACCTGTATTTCGCTGTGGCTTCAACATCAACGTTCGGGACTTTGAGGGGAAGGTTGTGAAGACTTCGGTTGTGTTCCACCAACTGGGGACAGCCATGCCTATGTCGGTGGAGGAAGGGCCTGAGTGCCAGGGACCTGTGGTTGACAGGCGCTGCCCTCGATGTGGTCATGAAGGAATGGCATACCACACCAGACAGATGCGTTTCAGCCGATGAAGGGCAAACCTGTCTTCTACACCTGTACCAACTGCAAGTTCAGGAGAAGGAAGACTCT ACGCGT ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC235838 representing NM_001278785 Red=Cloning site Green=Tags(s) MSVMDLANTCSSFQSDLDFCSDCGSVLPLPGAQDTVTCIRCGFNINVRDFEGKVVKTSVVFHQLGTAMPMSVEEGPECQGPVVDVRRPCRCHEGMAYHTRQMRSADEGQTVFYTCNCKFQEKEDS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001278785



[View online »](#)

ORF Size:	378 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001278785.1</u>
RefSeq Size:	760 bp
RefSeq ORF:	381 bp
Locus ID:	30834
UniProt ID:	<u>Q9P1U0</u>
Cytogenetics:	6p22.1
MW:	13.9 kDa