

Product datasheet for RC235843

TROAP (NM_001278324) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TROAP (NM_001278324) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TROAP
Synonyms: TASTIN
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC235843 representing NM_001278324
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACCACCCGGCAAGCCACGAAGGATCCCTCCTCCGGGGTGTATCTCCTACCCCTAGCAAGATTCCGG
 TACGCTCTCAGAAACGCACGCCTTTCCCACTGTTACATCGTGCCTGGACAGGAGAACCAAGATCC
 AAGGAGATGGGTGCAGAAACCACCGCTCAATATTCAACGCCCCCTCGTTGATTACGAGGCCCCAGGCCG
 AAAGCCAGGCACCAGGCAGAGACATCACAAGATTGGTGGGGATCAGTCAGCCTCGGAACCCCTTGAAG
 AGCTCAGGCCTAGCCCTAGGGGTCAAAATGTGGGGCTGGGCCCTGCCAGACAGGTGCCATGGCTTG
 TACCTGTAATCACAGCTACTCAGGAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC235843 representing NM_001278324
 Red=Cloning site Green=Tags(s)

MTTRQATKDLLRGVSPTPSKIPVRSQKRTPFPTVTSCAVDQENQDPRRWQKPLNIQRPLVDSAGPRP
 KARHQAETSQRLVGISQPRNPLEELRPSRPGQNVGPPAQTGAMACTCNHSSGG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_001278324

ORF Size: 378 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001278324.2](#)

RefSeq Size: 678 bp

RefSeq ORF: 381 bp
 Locus ID: 10024
 UniProt ID: Q12815
 Cytogenetics: 12q13.12
 MW: 14.2 kDa

Gene Summary: Could be involved with bystin and trophinin in a cell adhesion molecule complex that mediates an initial attachment of the blastocyst to uterine epithelial cells at the time of the embryo implantation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC235843