

## Product datasheet for **RC235861A1V**

### Human C3orf14 (NM\_001291942) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human C3orf14 (NM\_001291942) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** C3orf14  
**Synonyms:** HT021  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC235861 representing NM\_001291942  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGACTTCCTTGTGCTCAAGAAATTCGCCTTCTAAAAGACATGAAGAAATAGTATCACAAAGATTAA  
TGTTACTTCAACAAATGGAGAATAAATTGGGTGATCAACACACAGAAAAGGCATCTCAACTCCAACTGT  
TGAGACTGCTTTAAAAGGAACCTTAGTCTTTAAAGGATATAGAAGCAGCAGAAAAGTCACTACAGACC  
AGGATTCACCCACTTCCACGGCCTGAGGTGTTTCTCTTGAGACTCGTTACTGGGCATCAGTAGAAGAAT  
ATATCCCAAATGGGAACAGTTTCTTTAGGAAGAGCACCATATCCTTTTGCTGTTGAAAATCAAATGA  
AGCAGAAAATACCATTCAAATGAGGCACAGCGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC235861 representing NM\_001291942  
Red=Cloning site Green=Tags(s)

MTSLFAQEIRLSKRHEEIVSQRLMLLQQMENKLGDQHTEKASQLQTVETAFKRNLSELLKDIEAAEKSLQT  
RIHPLRPEVVSLETRYWASVEEYIPKWEQFLLGRAPYPFAVENQNEAENTIQNEAQR

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001291942



[View online »](#)

<b>ORF Size:</b>	384 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001291942.1</a></u>
<b>RefSeq Size:</b>	3320 bp
<b>RefSeq ORF:</b>	387 bp
<b>Locus ID:</b>	57415
<b>UniProt ID:</b>	<u><a href="#">Q9HBI5</a></u>
<b>Cytogenetics:</b>	3p14.2
<b>MW:</b>	15 kDa