

Product datasheet for **RC235897**

OPALIN (NM_001284324) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OPALIN (NM_001284324) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: OPALIN
Synonyms: HTMP10; TMEM10; TMP10
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235897 representing NM_001284324
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCGGATGACGTCTCTCTGTACAGGTGGAAAGAAACGGACTGTGGGCCCTCTTTGGATTAGCGG
CGGGCATAACCATTGCTGGTGGCCACAGCCCTGCTGGTGGCTTTACTATTTACTTTGATTCACCGAAGAAG
AAGCAGCATTGAGCCATGGAGGAAAGTGACAGACCATGTGAAATTTAGAAATTGATGACAATCCCAAG
ATATCTGAGAATCCTAGGAGATCACCCACACATGAGAAGAATACGATGGGAGCACAAAGAGGCCACATAT
ATGTGAAGACTGTAGCAGGAAGCGAGGAACCTGTGCATGACCGTTACCGTCTACTATAGAAATGGAAAG
AAGGAGGGGATTGTGGTGGCTTGTGCCAGACTGAGCCTGGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235897 representing NM_001284324
Red=Cloning site **Green**=Tags(s)
MRMTSSPVTGGKETDCGPSLGLAAGIPLLVAALLVALLFTLIHRRRSSIEAMEESDRPCEISEIDNPK
ISENPRRSPTHEKNTMGAQEAHIYVKTVAGSEEPVHDIRYPTIEMERRRGLWWLVPRLSLE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001284324

ORF Size: 393 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001284324.1 , NP_001271253.1
RefSeq Size:	3345 bp
RefSeq ORF:	396 bp
Locus ID:	93377
UniProt ID:	Q96PE5
Cytogenetics:	10q24.1
Protein Families:	Transmembrane
MW:	15.1 kDa
Gene Summary:	Central nervous system-specific myelin protein that increase myelin genes expression during oligodendrocyte differentiation. Promotes oligodendrocyte terminal differentiation. [UniProtKB/Swiss-Prot Function]