

## Product datasheet for **RC235899**

### ZNF302 (NM\_001289190) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF302 (NM\_001289190) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** ZNF302  
**Synonyms:** HSD16; MST154; MSTP154; ZNF135L; ZNF140L; ZNF327  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC235899 representing NM\_001289190  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCTCAGGTGACATTTAGTGATGTGGCTATAGACTTCTCTCATGAAGAGTGGGCATGCCTAGATTCTG  
CTCAGAGGGACTTATACAAGGATGTGATGGTCCAGAATTATGAGAACCTGGTCTCTGTAGGTCTTCCGT  
AACTAAGCCATATGTGATCATGTTATTGGAGGATGGAAAAGGCCCTGGATGATGGAGAAAAAAGTCA  
AAAGCTTACCCATTTCTTTATCACACTCTGTTCTGCTTCTGTGAACCTTGGATTCTCTGCTCTATTTG  
AGCATTGTTCAAGTCACTGAAATATTTGAGTTGTCAGAACTATGTGTTTTCTGGGTGCTTCATTTCTT  
ATCCAATTCTCCTAATCCACTGTAGAAGCTTTTTTCAAGAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC235899 representing NM\_001289190  
Red=Cloning site Green=Tags(s)

MSQVTFSDVAIDFSHEEWACLDSAQRDL YKDV MVQNYENL VSVGLSVTKPYVIMLLEDGKEPMMMEKLS  
KAYPFPLSHSVPASVNFVGFSALEHCSEVTEIFELSELVFWLHFLSNPSTVEAFFKK

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI




[View online >](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001289190

ORF Size: 393 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001289190.2</a>
<b>RefSeq Size:</b>	2888 bp
<b>RefSeq ORF:</b>	396 bp
<b>Locus ID:</b>	55900
<b>UniProt ID:</b>	<a href="#">Q9NR11</a>
<b>Cytogenetics:</b>	19q13.11
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	15.4 kDa
<b>Gene Summary:</b>	This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2014]