

Product datasheet for **RC235909**

FAM107B (NM_001282700) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: FAM107B (NM_001282700) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: FAM107B
Synonyms: C10orf45; HITS
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC235909 representing NM_001282700
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCCGAGCCAGACTACATAGAAGATGACAATCCTGAACTCATTAGGCCTCAGAACTGATCAATCCTG
TAAAAACCTCCCGAACCATCAAGATCTTACAGAGAACTTCTTATGAATCAAAAAAGGGTCTTGCTCC
TCAGAACAACCAGAATTGCAGAAGGTGATGGAAAAAGAAAACGAGACCAAGTAATAAAGCAGAAGGAA
GAAGAAGCACAGAAGAAGAAATCTGACTTGAAAATAGAGCTATTAACCGCAGCAGAAAGTTGGAGCAGC
TTGAACTTGAGAAGCAGAAATTGCAAGAAGAGCAAGAAAATGCCCCGAGTTTGTGAAGGTGAAAGGCAA
TCTCAGGAGAACAGGCCAAGAAGTCGCCAAGCCAGGAGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235909 representing NM_001282700
Red=Cloning site Green=Tags(s)

MAEPDYIEDDNPQLIRPQKLINPVKTSRNHQDLHRELLMNQKRGLAPQNKPELQKVMKRKRQVQIKQKE
EEAQKKSDLEIELLKRQKLEQLELEKQKLQEEQENAPEFVKVKGNLRRRTGQEVAAQES

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001282700.2
RefSeq Size:	3257 bp
RefSeq ORF:	396 bp
Locus ID:	83641
UniProt ID:	Q9H098
Cytogenetics:	10p13
MW:	16 kDa