

Product datasheet for RC235927

EBPL (NM_001278636) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: EBPL (NM_001278636) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: EBPL
Synonyms: EBRP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC235927 representing NM_001278636
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGGCGCTGAGTGGGAGCTGGGGGCCGAGGCTGGCGTTTCGCTGCTGCTGTGCGCCGCGCTGCTGGCGG
 CGGGCTGCGCCCTGGGCCTGCGCCTGGGCCGCGGGCAGGGGGCGGCGGACCGGGGCGCTCATCTGGCT
 CTGCTACGACGCGCTGGTGCACCTCGCGCTGGAAGGCCCTTTGTCTACTTGTCTTTAGTAGGAAACGTT
 GCAAATCCGATGGCTTGATTGCTTCTTTATGAAAGAATATGGCAAAGCTGATGCAAGATGGGTTTATT
 TTGATCCAACCATTTGTGTCTGTGAAATTCTGACCGTCGCCCTGGATGGGTCTCTGGCATTGTTCCCTCAT
 TTATGCCATAGTCAAAGAAAAATATTACCGGGACAGAGTGGCCTGGGGTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235927 representing NM_001278636
 Red=Cloning site Green=Tags(s)

MGAEWELGAEAGGSLLLCAALLAAGCALGLRLGRQGGAADRGALIWLCYDALVHFALGPFVYLSLVGNV
 ANSDGLIASLWKEYGKADARWVYFDPTIVSVEILTVALDGSLALFLIYAIIVKEKYRDRVAWGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_001278636

ORF Size: 402 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001278636.1](#), [NP_001265565.1](#)

RefSeq Size: 1038 bp

RefSeq ORF: 405 bp
 Locus ID: 84650
 Cytogenetics: 13q14.2
 Protein Families: Transmembrane
 MW: 14.8 kDa

Gene Summary: Does not possess sterol isomerase activity and does not bind sigma ligands.
 [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC235927