

Product datasheet for **RC235928**

CCDC32 (NM_001289132) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CCDC32 (NM_001289132) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: CCDC32
Synonyms: C15orf57
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC235928 representing NM_001289132
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGAAATGTTTGAGAGCGCTGACTCTACAGCCACAAGATCTGGCCAGGATCTCTGGGCTGAAATTTGTT
 CCTGTCTGCCAAATCCTGAACAAGAAGATGGTGCCAACAATGCATTCTCAGACTCCTTTGTGGATTCTTG
 CCCTGAAGGTGAAGGCCAGAGGGAGGTGGCTGACTTTGCTGTCCAGCCAGCTGTAAAGCCTTGGGCTCCC
 TTGCAGGATTAGAAGTGTATTTAGCATCTCTAGAGAAGAAGCTAAGAAGAATCAAAGGTTTAAATCAGG
 AAGTGACTTCCAAGGACATGCTTCGAACCTCTGGCCCAAGCAAGAAGGAATGCTGGGATCGGTTCTCTCA
 GGAGAAGTTAGCTTCAGAGTTCTTTGTGGATGGACTTGATTCTGATGAGAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235928 representing NM_001289132
Red=Cloning site Green=Tags(s)


MKMFESADSTATRSGQDLWAEICSLPNPEQEDGANNAFSDSFVDSCEPGEQREVADFAVQPAVKPWAP
 LQDSEVYLASLEKKLRRIKGLNQEVTSKDLRLTLAQAKKECWDRFLQEKLA SEFFVDGLDSDER

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM 001289132

ORF Size: 402 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001289132.2](#)

RefSeq Size: 1642 bp

RefSeq ORF: 405 bp
 Locus ID: 90416
 UniProt ID: Q9BV29
 Cytogenetics: 15q15.1
 MW: 15.5 kDa

Product images:



Circular map for RC235928