

Product datasheet for **RC235959A1V**

Human Host cell factor C1 regulator 1 (HCFC1R1) (NM_001288665) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Host cell factor C1 regulator 1 (HCFC1R1) (NM_001288665) AAV Particle
Tag: Myc-DDK
Symbol: Host cell factor C1 regulator 1
Synonyms: HPIP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC235959 representing NM_001288665
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATCCTGCAGCAGCCCTTGCAGCGAGGCCCCAGGGAGGGGCCAGCGCCTCCCGGGGCCGCTTGG
GGGTGACTTGGGGCTGGACGCCAGCTCCCTCTCCGAGGAGCTGTGCCATGAGCACAAGCGGCGCCT
GGAGGAGGAGCAGGAGCCTCTGCGCAAGCAGTTCTGTCTGAGGAGAACATGGCCACCCACTTCTCTCAA
CTCAGCCTGCACAATGACCACCCTACTGCAGCCCCCATGACCTTCTCCCCAGCCCTGCCCCACTCA
GGAGCCCTTGCTCTGAGCTGCTTCTCTGGCGCTATCCTGGCAGCCTCATCCCTGAGGCCCTCCGTCTGCT
GAGGCTGGGGGACACCCCGAGTCCCCCTACCCTGCAACCCAGCTGGGGACATAATGGAGCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC235959 representing NM_001288665
Red=Cloning site Green=Tags(s)

MILQQPLQRGPQGAQRLPRAALGVTWGLDASSPLRGAVPMSTKRRLEEEQEPLRKQFLSEENMATHFSQ
LSLHNDHPYCSPMFTSPALPPLRSPCSELLLWRYPGSLIPEALRLLRLGDTSPPPYPATPAGDIMEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001288665



[View online »](#)

ORF Size:	414 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001288665.1</u> , <u>NP_001275594.1</u>
RefSeq Size:	944 bp
RefSeq ORF:	417 bp
Locus ID:	54985
UniProt ID:	<u>Q9NWW0</u>
Cytogenetics:	16p13.3
MW:	15.3 kDa