

Product datasheet for **RC236014**

ADAT2 (NM_001286259) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ADAT2 (NM_001286259) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ADAT2
Synonyms: DEADC1; dj20N2; dj20N2.1; TAD2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC236014 representing NM_001286259
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGTCTACAACAATGAAGTTGTAGGGAAGGGGAGAAATGAAGTTAACCAAACCAAAAATGCTACTCGAC
ATGCAGAAATGGTGGCCATCGATCAGGTCCTCGATTGGTGTCTGCAAAGTGGCAAGAGTCCCTCTGAAGT
ATTTGAACACACTGTGTTGTATGTCCTGTGGAGCCGTGCATTATGTGTGCAGCTGCTCTCCGCCTGATG
AAAATCCCGCTGGTTGTATATGGCTGTCAGAATGAACGATTTGGTGGTTGTGGCTCTGTTCTAAATATTG
CCTCTGCTGACCTACCAAACTGGGAGACCATTTTCAGTGTATCCCTGGATATCGGGCTGAGGAAGCAGT
GGAAATGTTAAAGACCTTCTACAAACAAGAAAATCCAAATGCACCAAAATCGAAAGTTCCGAAAAAGGAA
TGTCAGAAATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236014 representing NM_001286259
Red=Cloning site Green=Tags(s)
MVYNNEVVGKGRNEVNQTKNATRHAEMVAIDQVLDWCRQSGKSPSEVFEHTVLYVTVEPCIMCAAALRLM
KIPLVVYGCQNERFGGCGSVLNIASADLPNTGRPFQCIPGYRAEEAVEMLKTFYKQENPNAPKSKVRKKE
CQKS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001286259

ORF Size: 432 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001286259.1, NP_001273188.1</u>
RefSeq Size:	7368 bp
RefSeq ORF:	435 bp
Locus ID:	134637
UniProt ID:	<u>Q7Z6V5</u>
Cytogenetics:	6q24.2
MW:	16.5 kDa
Gene Summary:	Probably participates in deamination of adenosine-34 to inosine in many tRNAs. [UniProtKB/Swiss-Prot Function]