

## Product datasheet for **RC236022**

### **TMEM68 (NM\_001286661) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM68 (NM\_001286661) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TMEM68  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC236022 representing NM\_001286661  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAAAATACCAGAAGATGGACCAGCACTTATAATTTTTTATCATGGAGCTATTCCTATAGATTTTT  
ACTATTTTCATGGCTAAAAATTTTATACACAAAGGCAGAACTTGCCGAGTAGTAGCTGATCACTTTGTCTT  
TAAATTCAGGGTTTAGTTTACTGGATGTGTTTTGTGCTCTACATGGACCAAGAGAAAAATGTGT  
GAAATTCGAGGAGTGGCCACTTGTTAGCTATCTCACCAGGTGGAGTTCGAGAAGCCCTAATTAGTGATG  
AACTTATAACATCGTATGGGGTCATCGCAGAGGCTTTGCTCAGGTTGCAATTGATGCAAAGTGCCCAT  
TATTCCTATGTTTACACAAAATATTCGAGAAGGATTTAGATCACTTGGAGGAACAAACGAAGAATGCTGT  
TCAAGCTTTGAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC236022 representing NM\_001286661  
**Red=Cloning site Green=Tags(s)**

MEKIPEDGPALIIIFYHGAIPIDFYFMAKIFIHKGRTRVVDHDFVKIPGFSLLLDVFCALHGPREKCV  
EILRSGHLLAISPGGVREALISDETYNIVGHRRGFAQVAIDAKVPIIPMFTQNIREFRSLGGTNEECC  
SSFD

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001286661

**ORF Size:** 432 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


**RefSeq:** [NM\\_001286661.2](#)

**RefSeq Size:** 2037 bp

**RefSeq ORF:** 435 bp

Locus ID: 137695  
 Cytogenetics: 8q12.1  
 Protein Families: Transmembrane  
 MW: 16.6 kDa

**Product images:**



Circular map for RC236022