

Product datasheet for **RC236027**

C5orf66 (NM_001277348) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C5orf66 (NM_001277348) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C5orf66
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC236027 representing NM_001277348
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGATAACACCCCTGGCCCAAGGGTCACACAGAACCAAAAGGAAGCTACAGCCTAAGGAAGAAACCCCT
 GCTATGGGACTGTGGCTGCTCTGGACTGCGCAGGCTGCCAGGAAGGAATGGCGGAGGCAGCCTCCGTGGG
 CCTTGCACAGAGTGAGAGACCAACAATTATGGGCCATGTGTGCGCAGAGAGACACCCCATGACCGATAC
 GTAACAAGCATGTGCCCTACCTGCCGTGCCGATGGATGGAAATGCAATGGATTTCAGTCTCTTATCAT
 CAGCCAGGGCCAAGCAGTTTTTCACTGTCTTTCCAGAAGTTCTTCACACTTGTCTGCACCCAACTGG
 ACTATTACAGTGGATCACAACTTGGCAGGCAGGCCATGGCCTGCGCTGACAGCAGCTAGTGGTGGGGG
 ATGAAACAGAAACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236027 representing NM_001277348
 Red=Cloning site Green=Tags(s)

MGITPWPKGHTPEKGSYSLRKKPCYGTVAALDCAGCQEGMAEASVGLAQSESTTIMGHVCAERHHTDRY
 VTSMPYLPMPDGNAMDFSFLSSARAKQFFTVFSRSSHLAPQTGLLQWITNLGRPWPALTAASGAG
 MKQKQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001277348

ORF Size: 435 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001277348.1](#), [NP_001264277.1](#)

RefSeq Size: 2490 bp


RefSeq ORF: 438 bp

Locus ID: 100996485

Cytogenetics: 5q31.1

MW: 16 kDa

Product images:



Circular map for RC236027