

Product datasheet for **RC236047A1V**

Human CARHSP1 (NM_001278262) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CARHSP1 (NM_001278262) AAV Particle
Tag: Myc-DDK
Symbol: CARHSP1
Synonyms: CRHSP-24; CRHSP24; CSDC1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC236047 representing NM_001278262
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCATCTGAGCCTCCCCACCACCACAGCCCCCACCATCAAGCTTCAGTCGGGCTGCTGGACACCC
 CTCGGAGCCGTGAGCGCTCACCATCCCCTCTGCGGGGCAACGTGGTCCCAAGCCCACTGCCACTCGCCG
 GACGAGGACCTTCTCGGCGACGGTGC GGCTTCACAGGGCCCCGTCTACAAAGGAGTCTGCAAATGCTTC
 TGCCCGTCCAAGGGCCATGGCTTCATTACCCAGCTGATGGCGGCCCGACATCTCTCGACATCTCTG
 ATGTGGAAGGGGAGTATGTCCCAGTGGAAAGCGACGAGGTCACCTATAAAATGTGCTCCATCCCACCCAA
 GAATGAGAAGCTGCAGGCCGTGGAGGTCGTCACCTACCTGGCACCAGGCACCAAGCATGAGACCTGG
 TCTGGACATGTCATCAGCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236047 representing NM_001278262
Red=Cloning site Green=Tags(s)

MSSEPPPPQPPTHQASVGLLDTPRSRERSPSPLRGNVVPSPLPTRRTRTF SATVRASQGPVYKGVCKCF
 CRSKGHGFITPADGGPDIFLHISDVEGEYVPVEGDEVYKMCSSIPKNEKLQAVEVVITHLAPGTHKETH
 SGHVISS

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



[View online »](#)

ACCN:	NM_001278262
ORF Size:	441 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001278262.1, NP_001265191.1</u>
RefSeq Size:	2949 bp
RefSeq ORF:	444 bp
Locus ID:	23589
UniProt ID:	<u>Q9Y2V2</u>
Cytogenetics:	16p13.2
MW:	15.9 kDa