

Product datasheet for **RC236075**

TMEM85 (EMC4) (NM_001286420) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM85 (EMC4) (NM_001286420) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TMEM85
Synonyms:	PIG17; TMEM85
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC236075 representing NM_001286420 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACGGCCAGGGGGCCTGGTGGCTAACCGAGGCCGGCGCTTCAAGTGGGCCATTGAGCTAAGCGGGC
CTGGAGGAGGCAGCAGGGGTGGAAGTGACCGGGCAGTGGCCAGGGAGACTCGCTCTACCCAGTCGGTTA
CTTGGACAAGCAAGTGCCTGATACCAGCGTGCAAGAGACAGACCGGATCCTGGTGGAGAAGCGCTGCTGG
GACATCGCCTTGGGTCCCTCAAACAGATTCCCATGAATCTCTTCATCATGTACATGGCAGGCAATACTA
TCTCCATCTTCCCTACTATGATGGTGTGTATGATGGCCTGGCGACCCATTCAGGCACTTATGGCCATTTT
AGCCAAGAATGGAGTTCAGTGGTGGAGGACTGCTTTTGTGAACATGAGAAAGCAGCGCCTGGTCCCTATG
TATTTGGGTCTTATTACATCCTTCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC236075 representing NM_001286420 Red=Cloning site Green=Tags(s)
-------------------	---

MTAQGGLVANRGRFRKWAIELSGPGGSRGRSDRSGQGDSLVPVGYLDKQVPDTSVQETDRILVEKRCW
DIALGPLKQIPMNLFIYMAGNTISIFPTMMVCMMAWRPIQALMAISAKNGVQWWRTAFVNMKQRLVPM
YLGLIYILL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	SgfI-MluI
--------------------	-----------



[View online »](#)

Cloning Scheme:


ACCN: NM_001286420

ORF Size: 447 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001286420.2](#)


RefSeq Size: 904 bp

RefSeq ORF: 450 bp

Locus ID: 51234

UniProt ID: [Q5J8M3](#)
 Cytogenetics: 15q14
 Protein Families: Transmembrane
 MW: 17.2 kDa
 Gene Summary: May mediate anti-apoptotic activity.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC236075