

Product datasheet for **RC236076**

SLC22A18AS (NM_001302862) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SLC22A18AS (NM_001302862) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SLC22A18AS
Synonyms: BWR1B; BWSCR1B; ORCTL2S; p27-BWR1B; SLC22A1LS
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC236076 representing NM_001302862
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCGGTGTGCAGAGGGTGCTTGGTGGTTCTCTCCTGATGGCCCCGAGGGTCTGCAGCCTCCATCTGGC
CAGCAGAGGGCGCAGAAGGACTGCCTGGGCAGCTCGGACGTGACCGCCTGGAAGTGGTGTACAGCGTTCC
TGACAACGTTCCCGCCAAAACGGGTCCCGCCGCCCACTTGTGTGCAAGATAACTGGAAAATGTCTTTCT
GTGTGCTCCGAGGAGAATGCAAAGGCTGGTGGATGTAGTGCTTTCTCTACTGCTCTCTCAGCTGGGGG
CAAGAATGACAGGACGTGAACATGCACACAAGGGCCCGAACTCACGACCCCGACAGCGGTCTCCCCCG
CCCCCCAACCCGCGCTTGCAGGATTTAGGGCACTAGCACAGCACAGTCCACCCCTTGGGACTAGCACC
CCCTCCGCAGTCTGCTCTCTGCAGCAACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236076 representing NM_001302862
Red=Cloning site Green=Tags(s)

MRCAEGAWWFSPDGPAGSAASIWPAEAGAEGLPGQLGRDRLEVYVSPDNVPGQNGSRRPLVCKITGKCLS
 VCSEENAKAGGCSAFPLLLSQLGARMTGREHAHKPELTTPDSGLPRPPNPALAGFRALAQHSPPLGTST
 PSAVLLSAAT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001302862

ORF Size: 450 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001302862.2](#)

RefSeq Size: 1916 bp

RefSeq ORF: 453 bp


Locus ID: 5003

UniProt ID: Q8N1D0

Cytogenetics: 11p15.4

MW: 15.9 kDa

Product images:



Circular map for RC236076