

Product datasheet for **RC236086**

ZNF550 (NM_001277093) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF550 (NM_001277093) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ZNF550
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC236086 representing NM_001277093
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCGGAGACGAAGGACGCAGCGCAGATGTTGGTGACCTTCAAGGATGTGGCTGTGACCTTACCCGGG
AGGAGTGGAGACAGCTGGACCTGGCCAGAGGACCTGTACCGAGAGGTGATGCTGGAGACCTGTGGGCT
TCTGGTTTCACTAGGGCATCGGGTCCCAAACCAGAGTTGGTCCACCTGCTAGAGCATGGCAGGAGCTG
TGGATAGTGAAGAGAGGCCTCTACATGCTACCTGTGCAGGTGACAGAGCACAAGTTCATACCAGAGAGC
CAACCACTTACCCGCCAGTCTTATCTGAGCGGGCTTTCTCCGGGAAGCCTGACTCTGGAATCTTCAAC
CTCAAGTGATTGAGGTTGGGAGAGCTAGGGATGAGGAAGGGCTTTTGGAAATGCAGAAAGAAAACTG
CTTAAGGTTGAAAAGAGAAATTTAAGAAATTGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236086 representing NM_001277093
Red=Cloning site Green=Tags(s)

MAETKDAAQMLVTFKDVAVTFTREEWRQLDLAQRLLYREVMLETCGLLVSLGHRVPKPELVHLLHGHQEL
WIVKRGLSHATCAGDRAQVHTREPTTYPVLSERAFLRGLTLESSTSSDSRLGRARDEGLLEMQKEKL
LKVEKRNLNRC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001277093

ORF Size: 453 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001277093.2
RefSeq Size:	2015 bp
RefSeq ORF:	456 bp
Locus ID:	162972
UniProt ID:	Q7Z398
Cytogenetics:	19q13.43
Protein Families:	Transcription Factors
MW:	17.7 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]