

Product datasheet for RC236173

KRTAP9-6 (NM_001277331) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: KRTAP9-6 (NM_001277331) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: KRTAP9-6
Synonyms: KAP9.6; KRTAP9.6; KRTAP9L2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC236173 representing NM_001277331
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACCCACTGTTGCTCCCTGGCTGTCAGCCTACCTGCTGCAGGACCACTTGCTGCAGGACTACCTGCT
 GGCAGCCCACCATTGTGACCACCTGCAGCAGCACACCCTGCTGCCAGCCCTCCTGCTGTGTGTCCAGCTG
 CTGCCAGCCTTACTGCCACCCAATTGCTGTCAAACACCTGCTGCAGGACCACCTGCTGCCAGCCCACC
 TGTGTGACCAGCTGCTGCCAGCCTTCTGCTGCAGCACACCCTGCTACCAGCCCATCTGCTGTGGGTCCA
 GCTGCTGTGGCCAAACCAGCTGTGGGTCCAGCTGTGGCCAGAGCAGCTCCTGTGCACCTGTGTACTGCAG
 AAGAACCTGCTACCACCCACGACTGTCTGCCTGCCTGGTTGCCTCAACCAGAGCTGTGGATCCAGCTGC
 TGCCAGCCCTGCTACTGCCAGCCTGCTGTGTGTCCAGCTGCTGCCAGCATTCTTGTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236173 representing NM_001277331
 Red=Cloning site Green=Tags(s)


MTHCCSPGCQPTCCRTTCCRTTCWQPTIVTTCSSSTPCCQPSCCVSSCCQPYCHPTCCQNTCCRTTCCQPT
 CVTSCCQPSCCSTPCYQPICGSSSCGQTSCGSSCGQSSSCAPVYCRRTCYPHTTVCLPGLNQSCGSSC
 CQPCYCPACCVSSCCQHSCC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001277331

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001277331.1](#), [NP_001264260.1](#)

RefSeq Size: 483 bp

RefSeq ORF: 483 bp

Locus ID: 100507608


UniProt ID: A8MVA2

Cytogenetics: 17q21.2

MW: 17.3 kDa

Gene Summary: In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC236173