

Product datasheet for RC236178

ERI3 (NM_001301701) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ERI3 (NM_001301701) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ERI3
Synonyms:	PINT1; PRNPIP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC236178 representing NM_001301701 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGAGATTGAGTCTACCTTTCACATGTATGTCCAGCCTGTAGTCCATCCACAGCTTACCCATTCTGTA
 CAGAGCTACCGGGATTATTCAGCCATGGTGGATGGTCAGCCAAGCCTGCAGCAAGTGTGGAGAGGGT
 CGATGAATGGATGGCGAAGGAAGGCCTCTAGATCCAAACGTCAAGTCAATTTTGTACCTGTGGAGAC
 TGGGACTTAAAGTCATGCTCCAGGCCAGTCCAGTACTTGGGCTTGCCAGTGGCGGATTACTTCAAGC
 AGTGGATTAATCTGAAAAGGCTTACAGCTTCGCCATGGGCTGCTGGCCCAAGAATGGACTTCTAGACAT
 GAACAAGGGCCTCAGCCTGCAACACATAGGCCGGCCCCACAGCGCATTGACGACTGCAAGAATTGCC
 AACATCATGAAGACACTCGCCTATCGAGGCTTCATCTTCAAGCAGACATCGAAGCCGTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC236178 representing NM_001301701 Red=Cloning site Green=Tags(s)
-------------------	---

MEIESTFHMVQPVVHPQLTPFCTELTGIIQAMVDGQPSLQQVLERVDEWMAKEGLLDPNVKSIFVTCGD
 WDLKVMLPGQCQYGLPVADYFKQWINLKKAYSFAMGCWPKNGLLDNMKGLSLQHIGRPHSGIDDCNKIA
 NIMKTLAYRGFIFKQTSKPF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	SgfI-MluI
--------------------	-----------


[View online »](#)

Cloning Scheme:


ACCN: NM_001301701

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001301701.1](#), [NP_001288630.1](#)

RefSeq Size: 1135 bp

RefSeq ORF: 483 bp


Locus ID: 79033

UniProt ID: O43414

Cytogenetics: 1p34.1

MW: 18.6 kDa

Product images:



Circular map for RC236178