

## Product datasheet for **RC236247**

### **C1orf86 (FAAP20) (NM\_001282673) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C1orf86 (FAAP20) (NM\_001282673) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** FAAP20  
**Synonyms:** C1orf86; FP7162  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC236247 representing NM\_001282673  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

**ATG**GAGGCGGCGCGGAGGCCGGCTGGGGTTGAGCCGCCGAGGCCGCCCGCGGGCGGGCCTTCTG  
GCGGCCGCCCTGGTTTCTCTGGGGGTGATGAGCGGGAGCGGCTCTGGGCCGAGCTACTGCGCACGGT  
GAGCCCGAGCTGATCCTGGATCAGAGGTGCCTTCACTGCCCGCTTCCAGGACAGGAGCCAGGTGC  
GGCCCGAGCCCACTGAAGTCTTCACTGTCGGACCAAGACCTTTTCTGGACACCTTTCCGCCGACC  
TGTGGGGCCCGGCCGTTCTACCGGCTGCTTACGGGGCAGGAGGGCACCTGGAATCCCCGCCAGGTC  
CCTGCCCCAGCGCCCGCACCTGATCCCTGCAGGGCCCCAGGGTGGAGCAGCAGCCGTCTGTGGAGGGT  
GCCGCGCCCTGCGCAGCTGCCCATGTGCCAGAAGGAGTTCCGCCCCAGGTACTTCAACTCCCGCCCC  
ACCTGTGTCCCGCTGGCTCC

**ACGCGT**ACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC236247 representing NM\_001282673  
**Red=Cloning site Green=Tags(s)**

MEAARRPRLGLSRRRPPAGGPPSGRPFLLGGDERERLWAE LLRTVSP ELILDHEVPSLPAFPGQEP RC  
GPEPTEVFTVGP KTF SWTPFPDLWGPGRSYRLLHGAGGHLESPARSLPQR PAPDPCRAPRVEQQPSVEG  
AAALRSCPMCQKEFAPRYFN SRPHLCPAGS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001282673

ORF Size: 510 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                               |   |
|-------------------------------|---|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001282673.1</a></u> , <u><a href="#">NP_001269602.1</a></u>   |
| <b>RefSeq Size:</b>           | 1567 bp   |
| <b>RefSeq ORF:</b>            | 513 bp  |
| <b>Locus ID:</b>              | 199990  |
| <b>UniProt ID:</b>            | <u><a href="#">Q6NZ36</a></u>   |
| <b>Cytogenetics:</b>          | 1p36.33   |
| <b>MW:</b>                    | 19.1 kDa  |
| <b>Gene Summary:</b>          | Component of the Fanconi anemia (FA) complex required to recruit the FA complex to DNA interstrand cross-links (ICLs) and promote ICLs repair. Following DNA damage recognizes and binds 'Lys-63'-linked ubiquitin generated by RNF8 at ICLs and recruits other components of the FA complex. Promotes translesion synthesis via interaction with REV1.[UniProtKB/Swiss-Prot Function]  |