

Product datasheet for **RC236381**

ARL6 (NM_001278293) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ARL6 (NM_001278293) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ARL6
Synonyms: BBS3; RP55
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC236381 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGATTGCTAGACAGACTTTTCAGTCTTGCTTGGCCTGAAGAAGAAGGAGGTTTCATGTTTTGTGCCTTG
 GGCTAGATAATAGTGGCAAAACGACGATCATTAAACAACTTAAACCTTCAAATGCTCAATCTCAAATAT
 CCTTCCAACAATAGGATTCAGCATAGAGAAATCAAATCATCCAGTTTGCATTTACAGTGTGGACATG
 TCAGGTCAAGGAAGATACAGAAATCTCTGGGAACACTATTATAAAGAAGGCCAAGCTATTATTTTGTCA
 TTGATAGTAGTAGATTAAGAATGGTTGTGGCCAAAGAAGAACTCGATACTCTCTGAATCATCCAGA
 TATTAAACACCGTCGAATCCAATCTTATTCTTTGCAAATAAAATGGATCTTAGAGATGCAGTGACATCT
 GTAAAAGTGTCTCAGTTGCTGTGTTAGAGAACATCAAAGATAAACCCCTGGCATATTTGTCTAGTGATG
 CCATAAAGGAGAAGGCTTGCAAGAAGGTGTAGACTGGCTTCAAGATCAGATCCAGACTGTGAAGACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236381 protein sequence
 Red=Cloning site Green=Tags(s)

MGLLDRLSVLLGLKKKEVHVLCLGLDNSGKTTIINKLKPSNAQSQNILPTIGFSIEKFKSSSLSFTVFD
 SGQGRYRNLWEHYKQAIIFVIDSSDRLRMVVAKEELDTLLNHPDIKHRRIPILFFANKMDLRDAVTS
 VKVSQLCLENIKDKPWHICASDAIKGEGLQEGVDWLQDQIQTVKT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6318_f03.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001278293

ORF Size: 558 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278293.3](#)

RefSeq Size: 1531 bp

RefSeq ORF: 561 bp

Locus ID: 84100


UniProt ID: [Q9H0F7](#)

Cytogenetics: 3q11.2

MW: 21.1 kDa

Gene Summary: The protein encoded by this gene belongs to the ARF-like (ADP ribosylation factor-like) sub-family of the ARF family of GTP-binding proteins which are involved in regulation of intracellular traffic. Mutations in this gene are associated with Bardet-Biedl syndrome (BBS). A vision-specific transcript, encoding long isoform BBS3L, has been described (PMID: 20333246). [provided by RefSeq, Apr 2016]

Product images:



Circular map for RC236381