

## Product datasheet for **RC236510**

### PPIL6 (NM\_001286361) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PPIL6 (NM\_001286361) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** PPIL6  
**Synonyms:** bA425D10.6; dj919F19.1; PPIase; RSPH12  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC236510 representing NM\_001286361  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCAAGGCCGAGCCGTGCGGGCCCCGCACGCTAGGTGCGGCTCGCCGTGCTGCCGAGCGGCCGC  
TGCAGGTGAAGGTGGTGGGCTCTCAGCTGCCCAACTTCAGATTGCGAAGAGCGCCGCTGAGGAACT  
CAAAAATGAAACCTGGGAATATTCTTCTCTGTGATTTCTTTTGTAAATGGTCAGTTTCTGGGTGATGCA  
TTGGATCTGCAGAAATGGGCCACGAGGTGTGGGATATAGTTGACATTAACCCTCTGCACCTTTATGACG  
CACTCACTGAGGATTTTTCCGCTAAGTTCTTAAGAGACACCAAGCATGATTTTCGTGTTTTGGACATTTG  
TATTGATCTTCTCCAATTGGAAGATTGATTTTTGAGCTATACTGTGATGTGTGCCAAAAACATGTAAA  
AATTTTCAGGTCTGTGCACAGGAAAAGCAGGGTTTTCTCAACGTGGCATAAGACTACATTACAAAAATT  
CCATTTTTCATCGAATAGTACAGAATGGCTGGATACAAGGAGGGATATAGTGATGGAAAAGGAGATAA  
TGGAGAGTCGATTTATGGTCCAACATTTGAAGGCAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC236510 representing NM\_001286361  
Red=Cloning site Green=Tags(s)

MARPQCGPPHARCGSPSLPERPLQVKVVGLFSCPWFQIAKSAEELKNETWEYSSSVISFVNGQFLGDA  
LDLQKWAHEVWDIVDIKPSALYDALTEDFSQKFLRDTKHDFVFLDICIDSSPIGRILIFELYCDVCPKTCK  
NFQVLCCTGKAGFSQRGIRLHYKNSIFHRIVQNGWIQGGDIVYKGDNGESIYGPTFEGN


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001286361

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001286361.1](#), [NP\\_001273290.1](#)

**RefSeq Size:** 3913 bp

**RefSeq ORF:** 600 bp

**Locus ID:** 285755

**UniProt ID:** [Q8IXY8](#)

**Cytogenetics:** 6q21

**MW:** 22.7 kDa

**Gene Summary:** Probable inactive PPIase with no peptidyl-prolyl cis-trans isomerase activity.[UniProtKB/Swiss-Prot Function]