

Product datasheet for RC236559

E2F6 (NM_001278276) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: E2F6 (NM_001278276) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: E2F6
Synonyms: E2F-6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC236559 representing NM_001278276
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGATCTTGTCAGATCTGCTCCCGGGGTATTCTTGACTTAAACAAGGTTGCAACGAACTGGGAGTCC
 GAAAGCGGAGAGTGTATGACATACCAATGTCTTAGATGGAATCGACCTCGTTGAAAAGAAATCCAAGAA
 CCATATTAGATGGATAGGATCTGATCTTAGCAATTTTGGAGCAGTCCCCAACAAAAGAAGCTACAGGAG
 GAACCTTCTGACTTATCAGCAATGGAAGATGCTTTGGATGAGTTAATTAAGGATTGTGCTCAGCAGCTGT
 TTGAGTTAACAGATGACAAAGAAAATGAAAGACTAGCATATGTGACCTATCAAGACATTCATAGCATTCA
 GGCTTCCATGAACAGATCGTCATTGCAGTTAAAGCTCCAGCAGAAACCAGATTGGATGTTCCAGCTCCC
 AGAGAAGACTCTATCACAGTGCACATAAGGAGCACCAACGGACCTATCGATGTCTATTTGTGTGAAGTGG
 AGCAGGGTCAGACCAGTAACAAAAGTCTGAAGGTGTCGGGACCTCTTCATCTGAGAGCACTCATCCAGA
 AGGCCCTGAGGAAGAAGAAAATCCTCAGCAAAGTGAAGAATTGCTTGAAGTAAGCAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC236559 representing NM_001278276
 Red=Cloning site Green=Tags(s)

MDLVRSAPGGILDNLNKVATKLGVRKRRVYDITNVLDGIDLVEKSKNHIRWIGSDL SNFGAVPQQKLLQE
 ELSDSLAMEDALDEL IKDCAQQLFEL TDDKENERLAYVTYQDIHSIQAFHEQIIVIAVKAPAEATRLDVPAP
 REDSITVHIRSTNGPIDVYLCEVEQGQTSNKRSEGVGTSSSESTHPEGPEEEENPQQSEELLEVSN


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001278276

ORF Size: 618 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001278276.2](#)

RefSeq Size: 3379 bp

RefSeq ORF: 621 bp

Locus ID: 1876

UniProt ID: [O75461](#)

Cytogenetics: 2p25.1

Protein Families: Transcription Factors

MW: 23.4 kDa

Gene Summary: This gene encodes a member of a family of transcription factors that play a crucial role in the control of the cell cycle. The protein encoded by this gene lacks the transactivation and tumor suppressor protein association domains found in other family members, and contains a modular suppression domain that functions in the inhibition of transcription. It interacts in a complex with chromatin modifying factors. There are pseudogenes for this gene on chromosomes 22 and X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]