

Product datasheet for **RC236641**

Lysophospholipase 1 (LYPLA1) (NM_001279358) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lysophospholipase 1 (LYPLA1) (NM_001279358) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LYPLA1
Synonyms:	APT-1; APT1; hAPT1; LPL-I; LPL1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC236641 representing NM_001279358 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGCGGCAATAACATGTCAACCCCGCTGCCCGCCATCGTGCCCGCCCGGGAAGGCCACCGCTGCGG
TGATTTTCCTGCATGGATTGGGAGATACTGGGCACGGATGGGCAGAAGCCTTTCAGGTATCAGAAGTTC
ACATATCAAATATATCTGCCCGCATGCGCCTGTTAGGCCTGTTACATTAATATGAACGTGGCTATGCCT
TCATGGTTTGATATTATTGGGCTTTCACAGATTACAGGAGGATGAATCTGGGATTAACAGGCAGCAG
AAAATATAAAAGCTTTGATTGATCAAGAAGTGAAGAATGGCATTCTTCTAACAGAATTATTTGGGAGG
GTTTTCTCAGGGAGGAGCTTTATCTTTATATACTGCCCTTACCACACAGCAGAACTGGCAGGTGTCAC
GCACTCAGTTGCTGGCTTCCACTTCGGGCTTCTTTCCACAGTGCACGGGGATTGTGACCTTTGGTTC
CCCTGATGTTGGTTCTCTTACGGTGGAAAACTAAAAACATTGGTGAATCCAGCCAATGTGACCTTTAA
AACCTATGAAGGTATGATGCACAGTTCGTGTCAACAGGAAATGATGGATGCAAGCAATTCATTGATAAA
CTCCTACCTCCAATTGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001279358.2
RefSeq Size:	2554 bp
RefSeq ORF:	651 bp
Locus ID:	10434
UniProt ID:	O75608
Cytogenetics:	8q11.23
Protein Pathways:	Glycerophospholipid metabolism
MW:	23.7 kDa
Gene Summary:	This gene encodes a member of the alpha/beta hydrolase superfamily. The encoded protein functions as a homodimer, exhibiting both depalmitoylating as well as lysophospholipase activity, and may be involved in Ras localization and signaling. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 4, 6, and 7. [provided by RefSeq, Jul 2013]